
# Instructions For Use

## Blue Sky Bio, LLC Dental Implant Systems

<b>English</b>	Blue Sky Bio, LLC Dental Implant Systems
<b>Le français</b>	Systèmes d'implants dentaires Blue Sky Bio, LLC
<b>Deutsch</b>	Blue Sky Bio, LLC-Dentalimplantatsysteme
<b>Italiano</b>	Sistemi di impianto dentale Blue Sky Bio, LLC
<b>Español</b>	Sistemas de implantes dentales Blue Sky Bio, LLC
<b>Svenska</b>	Blue Sky Bio, LLC Dentalt Implantssystem
<b>Ελληνικά</b>	Συστήματα οδοντικών εμφυτευμάτων Blue Sky Bio, LLC

**Blue Sky Bio, LLC**  
800 Liberty Dr  
Libertyville, IL 60048  
USA

**mdi Europa GmbH**  
Langenhagener Straße 71  
D-30855 Langenhagen  
Germany

	Catalog number		Do not reuse
	Caution, consult accompanying documents		Caution federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner
	Manufactured for		Sterilized using Gamma Radiation
	Do not use if package is damaged		Batch code
	Authorized Representative in the European Union		Use by date
	Non-sterile		Medical device
	Unique device identifier		Date of manufacture
	Single sterile barrier system with protective packaging inside		Name and address of the implantation healthcare institution / provider
	Date of implantation		Patient name or Patient ID
	Information website for patients		Electronic Instructions For Use
	MR Conditional		

# English (en) (For US Only)

## Blue Sky Bio, LLC Dental Implant Systems

### Instructions For Use

#### Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

3mm diameter implants are indicated only for maxillary lateral incisors and mandibular incisors.

Narrow diameter implants are not indicated to be used unsplinted in the posterior areas of the mouth.

#### Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy (Chrome, Cobalt or Molybdenum allergy for abutments manufactured from such alloys), cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants.

Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

#### Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, osseointegration failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

3mm diameter implants are recommended for single-unit use in maxillary lateral and mandibular incisors only, with a maximum torque limit of 40 Ncm and are not indicated for posterior areas of the mouth.

#### MRI Safety Information

While these products have not undergone MRI safety testing, literature analysis indicates that scanning a Blue Sky Bio, LLC Implant System poses no significant risks under specific conditions and are considered as "MR Conditionally Safe" under the following scanning conditions:

- Static magnetic field strengths of  $\leq 3$  Tesla
- Magnetic field gradients up to 30T/M (3000G/cm).
- RF Excitation:  
Circularly Polarized (CP) for all landmarks above umbilicus - No excitation restrictions for all landmarks below the umbilicus
- Operating Mode: Normal Operating Mode in the allowed imaging zone
- Whole-body specific absorption rates (SAR) of 2W/kg
- Maximum Head SAR: Not evaluated for head landmark
- Scan Duration: No specific constraints due to implant heating

#### Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

**Procedural Precautions, Prosthetics**

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

**Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants**

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. The orientation of the indexing feature of the abutment should match the orientation of the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

Original Blue Sky Bio NP stock and custom abutments are recommended for molars and wide diameter implants for significantly added strength.

Blue Sky Bio custom millable abutments are intended to be only customized by computerized equipment.

The surfaces mating with the implant should not be modified by any means.

The design parameters for customized abutments fabricated from Titanium Alloy Blank Abutments are:

Minimum wall thickness – 0.4 mm

Minimum post height for single-unit restorations – 4.0mm

4.0 mm Minimum gingival height – 0.5 mm

Maximum gingival height – 6.7 mm

Maximum angulation – 30°



**Compatibility:**

BIO   Max NP	BIO   Max NP is surgically and prosthetically compatible with the Nobel Active NP implant system.
BIO   Max RP	BIO   Max RP is surgically and prosthetically compatible with the Nobel Active RP implant system.
BIO   Internal Hex	BIO   Internal Hex is surgically and prosthetically compatible with the Zimmer/ZimVie Tapered Screw Vent and Implant Direct Internal Hex Legacy3 implant systems.
BIO   Quattro	BIO   Quattro is surgically and prosthetically compatible with the Straumann Bone Level (BL) implant system.
BIO   One Stage	BIO   One Stage is surgically and prosthetically compatible with the Straumann Tissue Level implant system.
BIO   Conus	BIO   Conus is surgically and prosthetically compatible with the Astra OsseoSpeed TX implant system.

**Implant Card Instructions for Completion**

An implant card has been provided with the implant device. It is intended to be completed by the healthcare provider; then giving to the patient.

1. Name of the patient or patient ID – To be completed by healthcare institution/provider
2. Date of implantation - To be completed by healthcare institution/provider
3. Name and address of the healthcare institution and provider - To be completed by healthcare institution/provider
4. Barcode to patient implant information – Present on Implant Card

Locations of text to be filled on the Left. Example of a completed implant card on the right with written information in blue.

Symbolling located on the device specific side of the implant card can be found on the first page of the IFU.

### **Follow Up Care**

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

### **Instructions for Base Abutments and for Custom abutments from Titanium Alloy Preformed Blanks**

#### **Indications for use**

Blue Sky Bio CAD-CAM Abutments are intended to be used in conjunction with Blue Sky Bio endosseous dental implants in the maxillary or mandibular arch to provide support for single-unit or multi-unit prosthetic restorations.

All digitally designed abutments for use with Blue Sky Bio CAD-CAM Abutments are intended to be sent to a Blue Sky Bio validated milling center for manufacture.

#### **Description**

Blue Sky Bio CAD-CAM Abutments include Titanium Alloy Base Abutments, Cobalt Alloy Base Abutments, and Titanium Alloy Blank Abutments. Titanium Alloy Base Abutments and Cobalt Alloy Base Abutments are provided in straight and 15° angled designs. Titanium Alloy Base Abutments and Cobalt Alloy Base abutments can support a zirconia superstructure or a final restoration. Cobalt Alloy Base Abutments may be cast to a final one-piece abutment.

#### **Fabrication of Suprastructure for Titanium Alloy and Cobalt Alloy Base Abutments**

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the superstructure using a computer software (CAD) the according to the design parameters below and the interproximal space available.

The recommend material for milling (CAM) the superstructure is Zirconia Y-TZP, conforming to ISO 13356 Implants for surgery - Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP).

#### **Superstructure Bonding**

The superstructure may be bonded to the abutment by the milling center, or may be bonded by the clinician. The following instructions are to be followed in either case.

Inspect the abutment for any defects. If any visible defects are present, do not use the abutment.

Grit blasting of the base abutment is not required for proper bonding due to the retentive features on the base.

Cleaning the base and the superstructure. Confirm there is no debris, grease, or other foreign matter present on any area to be cemented.

Multilink Hybrid Abutment Cement (Ivoclar Vivadent AG) is recommended. Refer to the manufacturer's instructions for proper mixing.

Apply cement to the contact surfaces of the base abutment and press the superstructure onto the abutment following the anti-rotational features.

Verify that both components are seated completely.

Remove excess cement immediately and refer to the manufacturer's instructions for proper curing procedures.

Sterilize abutment according to the steam autoclaving instructions.

#### **Cobalt Alloy Base Abutments Workflow for Casting**

**Cobalt Alloy Base Abutments (Straight) can be used as a base for a cast-to straight abutments only.**

**Angled Cobalt Alloy Base Abutments are not indicated for cast-to final abutment.**

The final one-piece abutment is fabricated using standard lost wax casting technique.

#### **Blue Sky Bio Titanium Alloy Preformed Blank Abutments**

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the final abutment using a computer software (CAD) the according to the design parameters below and the interproximal space available. The abutment is machined (CAM) using standard milling techniques.

#### **Design Parameters**

The design parameters for Titanium Alloy Base Abutments, Cobalt Alloy Base Abutments (including cast-to abutments) and for custom abutments fabricated from Titanium Alloy Preformed Blank Abutments are:

Minimum wall thickness – 0.4 mm

Minimum post height for single-unit restorations – 4.0 mm;

Minimum gingival height – 0.5 mm

Maximum gingival height – 6.7 mm

#### **Maximum angulation**

Base Abutments (Straight) – superstructure maximum angle – 30°

Base Abutments (Angled 15°) – superstructure maximum angle – 15° (total of 30°)

Final custom milled abutments – 30°

Cobalt Alloy cast-to abutments - 0°

Maximum divergence or angular correction between two implants is 60°

#### **Warning**

It is contraindicated to fabricate a zirconia superstructure and to cast a one-piece abutment, or to mill an abutment that exceed the design parameters.

### **Single Use Products:**

Single Use products must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

Single-use products are marked on the label by the "Do not reuse" symbol and should be disposed of after surgical usage.

### **Cleaning and Sterilization**

All implants are supplied sterile, and are for single use only prior to the labeled expiration date. Do not use implants if the packaging has been damaged or previously opened.

Abutments and screws are non-sterile and must be sterilized prior to use. See steam sterilization procedure below.

#### **Moist Heat (Steam) Sterilization Recommended Procedure:**

1. Place the device in a medical-grade paper/plastic steam-sterilization peel pouch compliant with EN ISO 11607-1 (or FDA-cleared for steam sterilization), sized to allow adequate seal margins and compatible with 121 °C gravity-displacement cycles.
2. Load Sterilizer (Gravity-Displacement) with pouches in a way that steam can circulate freely
3. Set sterilization parameters: 30-minute exposure at minimum of 121 °C (250 °F)
4. Run Sterilization cycle for the full 30-minute exposure at the set temperature
5. Allow 15 minutes of drying time before opening the sterilizer

Licensed Dental Clinicians are responsible for validating and routinely controlling their sterilization process and for ensuring their sterilization equipment is installed, maintained, and routinely monitored in accordance with the sterilizer manufacturer's instructions and applicable professional standards and local regulations (e.g., ANSI/AAMI ST79 and CDC guidance in the U.S.; EN ISO 17665 and national guidance in the EU).

#### **Method of Supply**

Blue Sky Bio, LLC implants are made out of medical grade titanium alloy (ASTM F136).

Blue Sky Bio, LLC abutments are made out of medical grade titanium alloy (ASTM F136), chrome-cobalt alloy (ASTM F1537) or gold alloy.

#### **Caution**

The sale of this device is restricted to, or by the order of, licensed physicians or dentists.

#### **Disposal**

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

Please report any serious incident to the manufacturer

Blue Sky Bio, LLC  
800 Liberty Drive  
Libertyville, IL 60048  
USA


# Instructions For Use

Blue Sky Bio, LLC Dental Implant Systems  
(en)

  
Blue Sky Bio, LLC  
800 Liberty Dr  
Libertyville, IL 60048  
USA

  
mdi Europa GmbH  
Langenhagener Straße 71  
D-30855 Langenhagen  
Germany

0297

	Catalog number		Do not reuse
	Caution, consult accompanying documents		Caution federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner
	Manufacturer		Sterilized using Gamma Radiation
	Do not use if package is damaged		Batch code
	Authorized Representative in the European Union		Use by date
	Non-sterile		Medical device
	Unique device identifier		Date of manufacture
	Single sterile barrier system with protective packaging inside		Name and address of the implantation healthcare institution / provider
	Date of implantation		Patient name or Patient ID
	Information website for patients		Electronic Instructions For Use
	MR Conditional		

# English (en) (For UK, Europe, Australia, and Canada)

## Blue Sky Bio, LLC Dental Implant Systems

### Instructions For Use

#### Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

3mm diameter implants are indicated only for maxillary lateral incisors and mandibular incisors.

Narrow diameter implants are not indicated to be used unsplinted in the posterior areas of the mouth.

#### Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy, cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants.

Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

#### Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, osseointegration failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

3mm diameter implants are recommended for single-unit use in maxillary lateral and mandibular incisors only, with a maximum torque limit of 40 Ncm and are not indicated for posterior areas of the mouth.

#### MRI Safety Information

While these products have not undergone MRI safety testing, literature analysis indicates that scanning a Blue Sky Bio, LLC Implant System poses no significant risks under specific conditions and are considered as "MR Conditionally Safe" under the following scanning conditions:

- Static magnetic field strengths of  $\leq 3$  Tesla
- Magnetic field gradients up to 30T/M (3000G/cm).
- RF Excitation:  
Circularly Polarized (CP) for all landmarks above umbilicus - No excitation restrictions for all landmarks below the umbilicus
- Operating Mode: Normal Operating Mode in the allowed imaging zone
- Whole-body specific absorption rates (SAR) of 2W/kg
- Maximum Head SAR: Not evaluated for head landmark
- Scan Duration: No specific constraints due to implant heating

#### Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

**Procedural Precautions, Prosthetics**

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

**Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants**

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. The orientation of the indexing feature of the abutment should match the orientation of the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

Original Blue Sky Bio NP stock and custom abutments are recommended for molars and wide diameter implants for significantly added strength.

Blue Sky Bio custom millable abutments are intended to be only customized by computerized equipment.

The surfaces mating with the implant should not be modified by any means.

The design parameters for customized abutments fabricated from Titanium Alloy Blank Abutments are:

Minimum wall thickness – 0.4 mm

Minimum post height for single-unit restorations – 4.0mm

4.0 mm Minimum gingival height – 0.5 mm

Maximum gingival height – 6.7 mm

Maximum angulation – 30°

**Compatibility:**

BIO   Max NP	BIO   Max NP is surgically and prosthetically compatible with the Nobel Active NP implant system.
BIO   Max RP	BIO   Max RP is surgically and prosthetically compatible with the Nobel Active RP implant system.
BIO   Internal Hex	BIO   Internal Hex is surgically and prosthetically compatible with the Zimmer/ZimVie Tapered Screw Vent and Implant Direct Internal Hex Legacy3 implant systems.
BIO   Quattro	BIO   Quattro is surgically and prosthetically compatible with the Straumann Bone Level (BL) implant system.
BIO   One Stage	BIO   One Stage is surgically and prosthetically compatible with the Straumann Tissue Level implant system.
BIO   Conus	BIO   Conus is surgically and prosthetically compatible with the Astra OsseoSpeed TX implant system.

**Implant Card Instructions for Completion**

An implant card has been provided with the implant device. It is intended to be completed by the healthcare provider; then giving to the patient.

1. Name of the patient or patient ID – To be completed by healthcare institution/provider
2. Date of implantation - To be completed by healthcare institution/provider
3. Name and address of the healthcare institution and provider - To be completed by healthcare institution/provider
4. Barcode to patient implant information – Present on Implant Card

<p><b>BlueSkyBio</b> International Implant Card</p> <p> 1 _____</p> <p> 2 _____</p> <p> 3 _____</p> <p> <a href="http://www.blueskybio.com/pages/patientimplantinfo">www.blueskybio.com/pages/patientimplantinfo</a> 4 </p>	<p><b>BlueSkyBio</b> International Implant Card</p> <p> <u>Jane Doe</u></p> <p> <u>03/15/2021</u></p> <p> <u>XYZ Surgical Center</u></p> <p><u>123 Healthcare Lane, Chicago IL</u></p> <p><u>Dr. A.B. Curtis</u></p> <p> <a href="http://www.blueskybio.com/pages/patientimplantinfo">www.blueskybio.com/pages/patientimplantinfo</a> </p>
---	--

Locations of text to be filled on the Left. Example of a completed implant card on the right with written information in blue. Symboling located on the device specific side of the implant card can be found on the first page of the IFU.

### Follow Up Care

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

### Single Use Products:

Single Use products must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

Single-use products are marked on the label by the "Do not reuse" symbol and should be disposed of after surgical usage.

Single-use products must not be reprocessed after surgical use. Re-use of single-use products creates a potential risk of contamination or device failure, which may lead to injury, illness or death of the patient.

### Cleaning and Sterilization:

All implants are supplied sterile, and are for single use only prior to the labeled expiration date Do not use implants if the packaging has been damaged or previously opened.

Abutments and screws are supplied non-sterile and must be sterilized prior to use. See steam sterilization procedure below.

### Moist Heat (Steam) Sterilization Recommended Procedure:

1. Place the device in a medical-grade paper/plastic steam-sterilization peel pouch compliant with EN ISO 11607-1 (or FDA-cleared for steam sterilization), sized to allow adequate seal margins and compatible with 121 °C gravity-displacement cycles.
2. Load Sterilizer (Gravity-Displacement) with pouches in a way that steam can circulate freely
3. Set sterilization parameters: 30-minute exposure at minimum of 121 °C (250 °F)
4. Run Sterilization cycle for the full 30-minute exposure at the set temperature
5. Allow 15 minutes of drying time before opening the sterilizer

Licensed Dental Clinicians are responsible for validating and routinely controlling their sterilization process and for ensuring their sterilization equipment is installed, maintained, and routinely monitored in accordance with the sterilizer manufacturer's instructions and applicable professional standards and local regulations (e.g., ANSI/AAMI ST79 and CDC guidance in the U.S.; EN ISO 17665 and national guidance in the EU).

### Method of Supply

Blue Sky Bio, LLC implants are made out of medical grade titanium alloy (ASTM F136).

Blue Sky Bio, LLC abutments are made out of medical grade titanium alloy (ASTM F136), or gold alloy.

### Caution

The sale of this device is restricted to, or by the order of, licensed physicians or dentists. (Prescription Use Only)

### Disposal:

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

### For customers in the European Union:

To request a copy of the Summary of Safety and Clinical Performance (SSCP) for implantable devices, please send an email specifying the article number(s) to [info@blueskybio.com](mailto:info@blueskybio.com)


Please report any serious incident to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Os implantes da Blue Sky Bio, LLC são realizados de liga de titânio de categoria médica (ASTM F136).
Os pilares da Blue Sky Bio, LLC são realizados de liga de titânio ou ouro de categoria médica (ASTM F136).