Instructions For Use

Blue Sky Bio, LLC Dental Implant Systems

English	Blue Sky Bio, LLC Dental Implant Systems
Le français	Systèmes d'implants dentaires Blue Sky Bio, LLC
Deutsch	Blue Sky Bio, LLC-Dentalimplantatsysteme
Italiano	Sistemi di impianto dentale Blue Sky Bio, LLC
Español	Sistemas de implantes dentales Blue Sky Bio, LLC
Svenska	Blue Sky Bio, LLC Dental Implant Systems
Türkçe	Blue Sky Bio, LLC Dental Implant Sistemleri
Ελληνικά	Συστήματα οδοντικών εμφυτευμάτων Blue Sky Bio, LLC
Português	Sistemas de implante dentário da Blue Sky Bio, LLC

Blue Sky Bio, LLC
800 Liberty Dr
Libertyville, IL 60048
USA

EC REP

mdi Europa GmbH
Langenhagener Straße 71
D-30855 Langenhagen
Germany

CE 0297

REF

CATALOG NUMBER

DO NOT REUSE

CAUTION, CONSULT ACCOMPANYING DOCUMENTS

Rx Only

CAUTION FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER

MANUFACTURED BY

STERILE R

STERILIZED USING GAMMA RADIATION

DO NOT USE IF PACKAGE IS DAMAGED

LOT

BATCH CODE

EC REP

AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

USE BY DATE

elfu Indicator

ELECTRONIC INSTRUCTIONS FOR USE

NON-STERILE

MD

MEDICAL DEVICE

UDI

UNIQUE DEVICE IDENTIFIER

DATE OF MANUFACTURE

SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE

English (For US Only)

Blue Sky Bio, LLC Dental Implant Systems

Instructions For Use

Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy (Chrome, Cobalt or Molybdenum allergy for abutments manufactured from such alloys), cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants. Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, implant failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

MRI Safety Information

Blue Sky Bio implants and abutments have not been evaluated for safety and compatibility in the Magnetic Resonance (MR) Environment. Blue Sky Bio implants and abutments have not been tested for heating, migration, or image artifact in the MR environment. The safety of Blue Sky Bio implants and abutments in the MR environment is unknown. Scanning a patient with Blue Sky Bio implants and abutments may result in patient injury.

Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

Procedural Precautions, Prosthetics

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. If the abutment is modified at chair side, it is advisable to make an impression and make a preliminary model to identify any undercuts prior to making the final impression. A carbide or diamond bur can be used with copious irrigation to remove the undercuts. The indexing feature of the abutment (e.g. octagon or trilobe) should match the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

Instructions for Base Abutments and for Custom abutments from Titanium Preformed Blanks

Indications for use

Blue Sky Bio CAD-CAM Abutments are intended to be used in conjunction with Blue Sky Bio endosseous dental implants in the maxillary or mandibular arch to provide support for single-unit or multi-unit prosthetic restorations.

All digitally designed abutments for use with Blue Sky Bio CAD-CAM Abutments are intended to be sent to a Blue Sky Bio validated milling center for manufacture.

Description

Blue Sky Bio CAD-CAM Abutments include Titanium Base Abutments, Cobalt Base Abutments, and Titanium Blank Abutments. Titanium Base Abutments and Cobalt Base Abutments are provided in straight and 15° angled designs. Titanium Base Abutments and Cobalt Base abutments can support a zirconia superstructure or a final restoration. Cobalt Base Abutments may be cast to a final one-piece abutment.

Fabrication of Suprastructure for Titanium and Cobalt Base Abutments

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the superstructure using a computer software (CAD) the according to the design parameters below and the interproximal space available.

The recommend material for milling (CAM) the superstructure is Zirconia Y-TZP, conforming to ISO 13356 Implants for surgery - Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP).

Superstructure Bonding

The superstructure may be bonded to the abutment by the milling center, or may be bonded by the clinician. The following instructions are to be followed in either case.

Inspect the abutment for any defects. If any visible defects are present, do not use the abutment.

Grit blasting of the base abutment is not required for proper bonding due to the retentive features on the base.

Cleaning the base and the superstructure. Confirm there is no debris, grease, or other foreign matter present on any area to be cemented.

Multilink Hybrid Abutment Cement (Ivoclar Vivadent AG) is recommended. Refer to the manufacturer's instructions for proper mixing.

Apply cement to the contact surfaces of the base abutment and press the superstructure onto the abutment following the anti-rotational features.

Verify that both components are seated completely.

Remove excess cement immediately and refer to the manufacturer's instructions for proper curing procedures.

Sterilize abutment according to the steam autoclaving instructions.

Cobalt Base Abutments Workflow for Casting

Cobalt Base Abutments (Straight) can be used as a base for a cast-to straight abutments only.

Angled Cobalt Base Abutments are not indicated for cast-to final abutment).

The final one-piece abutment is fabricated using standard lost wax casting technique.

Blue Sky Bio Titanium Preformed Blank Abutments

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the final abutment using a computer software (CAD) the according to the design parameters below and the interproximal space available. The abutment is machined (CAM) using standard milling techniques.

Design Parameters

The design parameters for Titanium Base Abutments, Cobalt Base Abutments (including cast-to abutments) and for custom abutments fabricated from Titanium Preformed Blank Abutments are:

Minimum wall thickness – 0.4 mm

Minimum post height for single-unit restorations – 4.0 mm;

Minimum gingival height – 0.5 mm

Maximum gingival height – 6.7 mm

Maximum angulation

Base Abutments (Straight) – superstructure maximum angle – 30°

Base Abutments (Angled 15°) – superstructure maximum angle – 15° (total of 30°)

Final custom milled abutments – 30°

Cobalt cast-to abutments - 0°

Maximum divergence or angular correction between two implants is 60°

Warning

It is contraindicated to fabricate a zirconia superstructure and to cast a one-piece abutment, or to mill an abutment that exceed the design parameters.

Follow Up Care

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

Sterility

All implants are supplied sterile, and are for single use only prior to the labeled expiration date (if applicable). Do not use implants if the packaging has been damaged or previously opened. Abutments and instruments are supplied non-sterile and must be cleaned and sterilized prior to use. Steam Sterilization procedure: Place the autoclave pouch containing item to be sterilized into the autoclave and follow the specific instructions provided by the manufacturer for pouched items. In general, a pouch must be sterilized by heating for 30 minutes at 250 deg F (121 deg C).

Cleaning instructions for reusable instruments

Rinse with cool-to-lukewarm water for two-and-one-half minutes. For all parts place in an ultrasonic cleaner with an enzymatic detergent diluted with tap water per the manufacturer's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes.

Method of Supply

Blue Sky Bio, LLC implants are made out of medical grade titanium alloy.

Blue Sky Bio, LLC abutments are made out of medical grade titanium alloy, chrome-cobalt alloy or gold alloy.

Caution

The sale of this device is restricted to, or by the order of, licensed physicians or dentists.

English (For UK, Europe, Australia, and Canada) Blue Sky Bio, LLC Dental Implant Systems

Instructions For Use

Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy, cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants. Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, implant failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

Procedural Precautions, Prosthetics

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. If the abutment is modified at chair side, it is advisable to make an impression and make a preliminary model to identify any undercuts prior to making the final impression. A carbide or diamond bur can be used with copious irrigation to remove the undercuts. The indexing feature of the abutment (e.g. octagon or trilobe) should match the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not

be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

Follow Up Care

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

Sterility

All implants are supplied sterile, and are for single use only prior to the labeled expiration date (if applicable). Do not use implants if the packaging has been damaged or previously opened. Abutments and instruments are supplied non-sterile and must be cleaned and sterilized prior to use. Steam Sterilization procedure: Place the autoclave pouch containing item to be sterilized into the autoclave and follow the specific instructions provided by the manufacturer for pouched items. In general, a pouch must be sterilized by heating for 30 minutes at 250 deg F (121 deg C).

Cleaning instructions for reusable instruments

Rinse with cool-to-lukewarm water for two-and-one-half minutes. For all parts place in an ultrasonic cleaner with an enzymatic detergent diluted with tap water per the manufacture's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes.

Method of Supply

Blue Sky Bio, LLC implants and abutments are made out of medical grade titanium or gold alloy.

Caution

The sale of this device is restricted to, or by the order of, licensed physicians or dentists.

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Português **Sistemas de implante dentário da Blue Sky Bio, LLC**
Instruções para Pilares de Restauração dos Implantes LLC da Blue Sky Bio
Os implantes e pilares da Blue Sky Bio, LLC são realizados de liga de titânio ou ouro de categoria médica.