Instructions For Use

Blue Sky Bio, LLC Single-Use Dental Implant Accessories

English	Blue Sky Bio, LLC Dental Implant Systems
Le français	Systèmes d'implants dentaires Blue Sky Bio, LLC
Deutsch	Blue Sky Bio, LLC-Dentalimplantatsysteme
Italiano	Sistemi di impianto dentale Blue Sky Bio, LLC
Español	Sistemas de implantes dentales Blue Sky Bio, LLC
Svenska	Blue Sky Bio, LLC Dental Implantssystem
Ελληνικά	Συστήματα οδοντικών εμφυτευμάτων Blue Sky Bio, LLC
Português	Sistemas de implante dentário da Blue Sky Bio, LLC

Blue Sky Bio, LLC
800 Liberty Dr
Libertyville, IL 60048
USA

mdi Europa GmbH
Langenhagener Straße 71
D-30855 Langenhagen
Germany

	Catalog number		Do not reuse
	Caution, consult accompanying documents		Caution federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner
	Manufactured for		
	Do not use if package is damaged		Batch code
	Authorized Representative in the European Community		
	Non-sterile		Medical device
	Unique device identifier		Date of manufacture
	Information website for patients		Electronic Instructions For Use

Indications For Use

Blue Sky Bio single-use restorative components are intended for use by licensed dental professionals in clinical and/or laboratory settings to:

- Transfer the three-dimensional position and orientation of an endosseous dental implant (or abutment) to an impression or working cast.
- Assist in the fabrication, fitting, and verification of prosthetic restorations supported by Blue Sky Bio implants.
- Assist in the fabrication of surgical drill guides for patients who have reached skeletal maturity

Device-specific intended uses

- **Analogs (implant/abutment replicas):** Used in impressions or casts to represent the implant/abutment position and facilitate prosthodontic laboratory procedures.
- **Impression transfers:** A device that duplicates the shape and position of the abutment interface of a dental implant in the impression to accurately position the implant analog
- **Guide tubes:** Assist with inserting drills during osteotomy as part of guided surgery workflows.
- **Waxing Sleeve, Space and Casting Sleeves:** For laboratory to facilitate direct waxing in preparation of creating a casting for a restoration or an abutment
- **Snap Inserts and Block-Out Spacers:** Use to replicate the position of overdenture-abutments in dentures.

Instructions

Impression transfers

1. Select correct platform and transfer type (open/closed-tray).
2. Seat and verify engagement and tighten fixation screw.
3. Make impression per standard clinical technique.
4. For closed-tray: remove transfer from implant, and reseat into impression.
5. For open tray: unscrew transfer before removal of impression tray.
6. Submit to dental laboratory specifying implant system/platform and restoration plan.

Analogs

Mate with impression transfer in the intended orientation; ensure full seating and secure fixation before pouring/processing.

Guide tubes

Use in accordance with the guided surgery plan to guide drills during osteotomy.

Snap Inserts and Block Out Spacers

Mate with overdenture-abutments to replicate the position of overdenture-abutments in dentures.

Waxing/burnout sleeves ,space & casting sleeves

Shape the sleeves and add wax to obtain the desired geometry. Ensure that the sleeves are heated to a temperature that is adequate for complete burn-out without any remaining residues and fill the resulting cavity completely with an appropriate molten dental metal.

Contraindications

Use is contraindicated when any general contraindications to impression-taking or restorative procedures are present (e.g., uncontrolled oral infection; inadequate stability of dental implants; known material hypersensitivity to device constituents). Do not use with incompatible implant platforms.

Warnings

- **Single-use:** Single-Use products must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross-contamination. Single-use components are marked on the label by the “Do not reuse” symbol  and should be disposed of after use.
- **Small parts / aspiration hazard:** Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry.
- **Verification:** Always verify full seating and engagement of transfers/analogues; take a radiograph where clinically indicated.
- **Compatibility:** Use only with the indicated Blue Sky Bio platform and corresponding components.

Table 1. Platform Compatibility—Single-Use Restorative Components

Blue Sky Bio Family	Single-use components	Compatible platform
BIO Max	Impression transfers, analogs, sleeves	NobelActive® NP/RP
BIO One Stage	Impression transfers, analogs, sleeves	Straumann-ITI® Tissue Level
BIO Quattro	Impression transfers, analogs	Straumann-ITI® Bone Level
BIO Internal Hex	Impression transfers, analogs	Zimvie (Previously Zimmer TSV®) / Implant Direct Legacy®
BIO Conus 12	Impression transfers, analogs	Astra Tech OsseoSpeed™ TX

Precautions

- Inspect each component prior to use; do not use if damaged or contaminated.
- Components that have a fixation screw should be tightened with light finger pressure.

Cleaning and Sterilization

Cleaning prior to use is not required for any single-use restorative component.

Do not sterilize prior to use: Analogs, plastic impression transfers, waxing sleeves, plastic snap inserts, space sleeves, casting sleeves, and block-out spacer. Use as supplied.

Sterilize prior to use: Guide tubes and Metal impression transfers must be sterilized prior to use. See steam sterilization procedure below.

Moist Heat (Steam) Sterilization Recommended Procedure:

1. Place the device in a medical-grade paper/plastic steam-sterilization peel pouch compliant with EN ISO 11607-1 (or FDA-cleared for steam sterilization), sized to allow adequate seal margins and compatible with 121 °C gravity-displacement cycles.
2. Load Sterilizer (Gravity-Displacement) with pouches in a way that steam can circulate freely.
3. Set sterilization parameters: 30-minute exposure at minimum of 121 °C (250 °F).
4. Run Sterilization cycle for the full 30-minute exposure at the set temperature.
5. Allow 15 minutes of drying time before opening the sterilizer.

Licensed Dental Clinicians are responsible for ensuring their sterilization equipment is installed, maintained, and routinely monitored in accordance with the sterilizer manufacturer's instructions and applicable local regulations (e.g., ANSI/AAMI ST79 and CDC guidance in the U.S.; EN ISO 17665 and national guidance in the EU).

Method of Supply

Blue Sky Bio, LLC single-use restorative components are supplied non-sterile unless stated otherwise on labeling. Components are made from metals (e.g., titanium alloy, stainless steel or aluminum) for metal transfers/sleeves; polymeric resins for plastic transfers/sleeves and burnout components.

Caution

Federal law restricts this device to sale by or on the order of a licensed Clinician (Rx Only).

Disposal

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

The hard copy of the IFU can be made available immediately through email, fax or direct mail free of cost.

Complaints and Adverse Events & Reporting

Please report any complaints or serious incidents to Blue Sky Bio customer service at (718) 376-0422 or <https://blueskybio.com/pages/support>.

Please report any serious incident to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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