dental.formlabs.com

formlabs 😽 | dental

Dentistry Made Easier



Materials and Applications

Unlock Every Dental Application With the Widest Range of Materials in the Industry

Formlabs offers the widest range of resins tailored to meet the unique needs of dental users, so you can produce safe and effective end-use appliances and restorations that can be patient-matched and produced at scale.

All Formlabs Dental resins are formulated and manufactured at our state-of-the-art facilities within a robust Quality Management System that is **ISO 13485 Certified, FDA-regulated, and EU MDR Compliant.** Furthermore, our biocompatible resins are manufactured to the highest global standards in our ISO Class 8 Clean Room.



Orthodontic Applications

⁶ Models for Diagnostics and Thermoforming Appliances



Models in Record Time for Diagnosis or Aligner Production

⁷ Indirect Bonding Trays

Accurate Bonding That's Easier to Plan and Use



⁸ Occlusal Splints and Guards

Appliances

For Long-Lasting Rigid and Flexible



Restorative Applications

¹⁰ Wax-Up Models

Your Digital Design With a Smoother Surface Finish



¹¹ Restorative Models

Soft and Rigid Components for Better Prosthetic Planning



¹² Temporary and Permanent Restorations

State-of-the-Art Ceramic-Filled Restorations



¹³ Direct Composite Restoration Guides

Perfectly Replicated Digital Designs With Faster Workflows



¹⁴ Custom Impression Trays Directly Printed for a Quicker Turnaround



¹⁵ Surgical Guides

For Highly Accurate and Predictable Implant Placement



¹⁶ Temporary All-on-X Appliances

~

For Reliable and Accurate Temporary Full-Arch Implant-Supported Restorations

¹⁷ Digital Dentures

Expanding Access to High-Quality Dentures in Custom Shades



¹⁸ Patterns for Casting and Pressing

Optimize Your Analog Workflow With Digital Accuracy



Formlabs Open Platform Solutions

Maximize the potential of your Formlabs 3D printer. Open Platform unlocks best-in-class third-party resins, custom print settings for tailored print performance, or the ability to experiment with any 405 nm photopolymer resin.

Orthodontic Applications



Models for Diagnostics and Thermoforming Appliances

FAST AND HIGHLY PRECISE MODELS IN RECORD TIME

Print 11 dental models in nine minutes, producing models for thermoforming clear aligners and retainers faster than ever. 3D print high-quality models in record time for case presentation and diagnosis. Streamline your workflow with PreForm Dental's Scan to Model feature for printable models in seconds.



Indirect Bonding Trays

ACCURATE BONDING: EFFORTLESSLY PLAN AND PRINT

Reduce chairtime and optimize bracket positioning by fabricating indirect bonding trays. IBT Flex Resin offers enhanced tear strength, translucency, and flexibility so you can give your patients a vastly improved experience.

	A CONTRACT OF THE OWNER OWNE
IBT Flex Resin	Fast Model Resin
Class I Biocompatible One-time use	Non-biocompatible
Indirect Bonding Trays Print high-quality indirect bonding trays with enhanced flexibility and tear strength	Models for Bracket Placement and Thermoformed Appliances Combine digital and analog procedures by printing models with brackets and thermoforming indirect bonding trays
100 microns	100 microns
Print time: 54 min	Print time: 25 min
Wash: 20 min	Wash: 5 min
Form Cure: 30 min at 70 °C Fast Cure: 5 min	Form Cure: 5 min Fast Cure: 1 min
\$/€ 2.50-4 per model	\$/€ 1-2 per model
~100	~100

Resin name

Classification

When to use

Suggested layer height

Workflow time For a full build platform

Cost per part

Parts per cartridge

Technical Data Sheet

Occlusal Splints and Guards

FOR LONG-LASTING RIGID OR FLEXIBLE APPLIANCES

Directly print rigid or flexible occlusal splints and guards with high accuracy, durability, optical transparency, discoloration resistance, and comfort for improved patient adoption and long-term outcomes.

Dental LT Clear Resin (V2)	Dental LT Comfort Resin
510(k) Cleared (US) / Class IIa (EU) Biocompatible Long-term use	510(k) Cleared (US) / Class IIa (EU) Biocompatible Long-term use
Rigid Occlusal Splints and Nightguards A stiff, highly durable, and fracture-resistant material that polishes to high optical transparency and resists discoloration over time	Flexible Occlusal Splints and Nightguards A flexible and durable material for comfortable long-term splints, and occlusal guards that are easily polished to high optical transparency
100 microns	100 microns
Print time: 43 min	Print time: 43 min
Wash: 15 min	Wash: 10 min
Form Cure: 60 min at 60 °C Fast Cure: 6 min	Form Cure: 20 min at 60 °C Fast Cure: 5 min
\$/€ 3.50-4.50 per part	\$/€ 4-5 per part
~100	~100

Resin name

Classification

When to use

Suggested layer height

Workflow time For a full build platform

Cost per part

Parts per cartridge

Technical Data Sheet

Form Fast C

~100



Restorative Applications



Wax-Up Models

YOUR DIGITAL DESIGN WITH A SMOOTHER SURFACE FINISH

From motivational mock-ups to testing aesthetics and function, users can accurately transfer digital smile designs with silicone indexes produced upon models 3D printed with Fast Model Resin, Precision Model Resin, or Grey Resin.



Restorative Models

SOFT AND RIGID COMPONENTS FOR BETTER PROSTHETIC PLANNING

Create high-accuracy models with flexible gingiva masks in combination with rigid dental models for more accurate planning of implant prosthetics. Customize your models with soft tissue printed in different shades.

	Store I	200 PG
	S A	
	Precision Model Resin	Soft Tissue Starter Pack
	Non-biocompatible	Non-biocompatible
	Light sand color with matte surface	Customizable dark, medium, and light pink shades
	Crisp and Accurate Models for Restorative Cases Consistently create restorative models with >99% of printed surface area within 100 µm of the digital model	Soft Tissue for Implant Models and Gingiva Masks Create flexible gingiva masks for use in combination with rigid dental models. Confidently check implant prosthetics by adding removable soft tissue components to your model production
height	50 microns	100 microns
arch printed form.	Print time: 1 h 25 min	Print time: ~1 h
	Wash: 5 min	Wash: 20 min
	Form Cure: 5 min at 35 °C Fast Cure: 30 sec	Form Cure: 10 min at 60 °C
	\$/€ 2-4 per model	\$/€ 1 per part
Sheet		

Resin name

Classification

Shades

When to use

Suggested layer height

Workflow time For one horseshoe arch printed flat to the build platform.

Cost per part

Technical Data Shee

Temporary and Permanent Restorations

STATE-OF-THE-ART CERAMIC-FILLED RESTORATIONS

Directly 3D print single-unit restorations and bridges with excellent marginal adaptation, strength, and aesthetics.



Resin name

Classification

VITA Shades

Compatibility

When to use

Suggested layer height

Workflow time For a single unit

Cost per part

Special hardware requirements

Technical Data Sheet

50 microns

Wash: 10 min

Form Cure: 30 min at 80 °C Fast Cure: 2x 4 min (flip halfway)

Special hardware requirements



Direct Composite Restoration Guides

PERFECTLY REPLICATED DIGITAL DESIGNS WITH FASTER WORKFLOWS

Equipped with high flexibility, tear resistance, and translucency for exceptional printing accuracy from a single restoration to a full smile. Save valuable technician time while delivering consistent, predictable outcomes.

	AL PARTY
Resin name	IBT Flex Resin
Classification	Class I (US) / Class I (EU) Biocompatible Short-term use
When to use	Direct Composite Restoration Guides With Enhanced Accuracy and Translucency From single units to full digital smile designs, 3D print flexible and tear-resistant translucent trays and guides that save you time and deliver consistent, predictable outcomes for both injectable and pressed composite techniques
Suggested layer height	50 microns
Workflow time For a single unit	Print time: ~2 h
	Wash: 20 min + 10 min soak or spray
	Form Cure: 30 min at 70 °C Fast Cure: 5 min
Cost per part	\$/€ 3-5 per guide
Technical Data Sheet	

Custom Impression Trays

DIRECTLY PRINTED FOR A QUICKER TURNAROUND

Directly print full impression trays for implants, dentures, crowns, bridges, and other comprehensive cases for reduced labor time, higher throughput, and consistent, accurate impressions with high-quality results.

Resin name	Custom Tray Resin
Classification	Class I (US) / Class I (EU) Biocompatible Short-term use
When to use	On-Demand Custom Trays For cases where you opt to use traditional PVS impressions, such as implant or fully edentulous cases, you can fabricate custom impression trays to assist in definitive impression taking
Suggested layer height	100 microns
Workflow time For a single unit	Print time: 39 min
	Wash: 10 min
	Form Cure: 30 min 60 ºC Fast Cure: 5 min
Cost per part	\$/€ ~6 per tray
Technical Data Sheet	

Surgical Guides

FOR HIGHLY ACCURATE AND PREDICTABLE IMPLANT PLACEMENT

Developed specifically for Formlabs printers and rigorously tested by dental specialists, our Surgical Guide Resin was designed to exceed demands in part quality, accuracy, and performance for better surgical outcomes.

Resin name	Surgical Guide Resin
Classification	Class I (US) / Class I (EU) Biocompatible Short-term use
When to use	Accurate and Autoclavable Surgical Guides You can fabricate precise and autoclavable surgical guides to aid in implant placement procedures of single or multi-unit implants
Suggested layer height	50 microns
Workflow Time For a single unit	Print time: 1 h 11 min
	Wash: 20 min
	Form Cure: 30 min 70 °C Fast Cure: 5 min
Cost per part	\$/€ 3-5 per guide
Technical Data Sheet	

15

Temporary All-on-X Appliances

TEMPORARY FULL-ARCH IMPLANT-SUPPORTED RESTORATIONS

Unlock 3D printed temporary full-arch implant-supported restorations in-house with our nano-ceramic filled Premium Teeth Resin. Equipped with optimal intraoral mechanical properties, fracture resistance, and accuracy, this material's aesthetics replicate a patient's natural smile, while saving valuable workflow time in post-processing and finishing steps.

Resin name	Premium Teeth Resin
Classification	Class II 510(k) Cleared (US) / Class IIa (EU) Biocompatible Temporary use (up to 12 months)
VITA Shades	BL A2 A3 B1
When to use	Accurate Temporary All-on-X Appliances 3D print temporary full-arch implant-supported restorations with optimal intraoral mechanical properties, fracture resistance, and accuracy
Suggested layer height	100 microns
Workflow time For a single unit	Print time: 26 min
	Wash: 10 min
	Form Cure: 30 min at 80 °C Fast Cure: 2x 4 min (flip halfway)
Cost per part	\$/€ 7-9 per part
Technical Data Sheet	

Digital Dentures

EXPANDING ACCESS TO HIGH-QUALITY DENTURES IN CUSTOM SHADES

Produce high-quality dentures in-house. Formlabs' Premium Teeth Resin and Denture Base Resin enable you to create full dentures in custom shades more easily and affordably than ever.



Resin name

Classification

VITA Shades

Compatibility

When to use

Suggested layer height

Workflow time For a single unit

Cost per part

Technical Data Sheet

Patterns for Casting and Pressing

FOR OPTIMIZING YOUR ANALOG WORKFLOW WITH DIGITAL ACCURACY

Castable Wax Resin provides casting and pressing patterns with sealed margins for accuracy, 20% wax for a clean burnout, and no curing requirements for a streamlined workflow.



Resin name

Classification

When to use

Your Next 3D Printing Idea:



Maximize Your Time, Money, and Throughput With Formlabs Dental 3D Printing Solutions

"We can actually do samevisit printing, where we can go from a scan to a printed model and fabricate an appliance very easily within 10 minutes so we have it before a patient leaves."

Dr. Christopher Baer DMD Baer Dental "It's a very easy set-up and very fast! From a user standpoint, meaning someone that is not trained to 3D print, this is a really easy set-up. It's so easy to use and intuitive"

Dr. Lisa Alvetro DDS Orthodontist & Owner Alvetro Orthodontics "Form 4B was a game changer for our entire production, because now we can produce parts two to five times faster than before."

Stephan Kreimer MDT Kreimer Dentallabor

