

Insert the Hex Driver within the 30Ncm torque ratchet

30 Ncm


30 Ncm


Insert driver into abutment and turn torque ratchet clockwise until torque of 30 Ncm is reached and the head of the ratchet releases Abutments in place

0


**Blue Sky Bio** 

Put silicon blockout ring on the abutment





# **Blue Sky Bio**



Remove enough material form the tissue side of the denture to ensure that it does not interfere with the housings.

Typically it is necessary to make the cavity in the denture 2-3mm wider and 1-2mm taller than the housing

**Blue Sky Bio** 

Insert denture and ensure that when it is fully seated it does not touch the housings

Insert resin of choice into the cavities in the denture

**Blue Sky Bio** 

Insert denture again over the housings and with the patient occluding on the denture

-



Let resin cure with denture in appropriate occlustion

**Blue Sky Bio** 

Void Coid Excess resin

> Remove denture from mouth and fill any voids with resin. Trim any excess resin around the housing and on the ridge area



Finished prosthesis with inserts in place

# **Blue Sky Bio**

Remove the insert you want to change with the sharp point of the removal instrument

**Blue Sky Bio** 

Place the insert with the desired retention into the denture with the opposite side of the removal instrument