




BlueSkyBio

Your Partner In Digital Dentistry

BIO | One Stage



2020

Product Catalog

BIO | One Stage

- Straumann Compatible

Valid January 2020

Table of contents	
Implants	2,3
Special Cover Screws	3
Solid Abutments and Restorative Kits	3-5
synOcta® Compatible Abutments	5
Gold Castable Abutments	5,6
Angled Abutments	6
Overdenture Abutments	6,7
Intraorally Preparable Abutments	7,8
Waxing Sleeves / Temporary Transfers	8
Abutment level analogs - Metal	8
Impression Transfers and Implant Level Analogs	9
Surgical Instruments	9,10
Safety Stop Surgical Kit	10,11
Other Surgical Instrumentation	12
Hammerless Crestal Sinus Lift Kit	12
BIO BlueTouch	12
Suture	13

Master Bone Expander Kit	13
Guided Sinus Lift Kit	13
Safety Lateral Sinus Kit	13
Guide Tubes	14-17
Terms & Conditions	18

Implants


Implants blasted with resorbable media and acid washed - Cover screw included.

Compatible with Straumann® ITI® implants, SynOcta® system, solid abutments, instrumentation and prosthetics.

Our new narrow platform uses the same prosthetics as regular platform.


For easy insertion, implants must be placed with a J-implant insertion driver and should not be used with other companies' mounts.

Narrow 3.3



Size	Part #
Ø3.3 x 10mm, 4.1mm Platform One Stage Implant	ION3310
Ø3.3 x 12mm, 4.1mm Platform One Stage Implant	ION3312

Narrow 4.1




Size	Part #
Ø4.1 x 10mm, 4.1mm Platform One Stage Implant	ION4110
Ø4.1 x 12mm, 4.1mm Platform One Stage Implant	ION4112

Regular 3.3



Size	Part #
Ø3.3 x 10mm, Regular Platform One Stage Implant	IOR3310
Ø3.3 x 12mm, Regular Platform One Stage Implant	IOR3312
Ø3.3 x 14mm, Regular Platform One Stage Implant	IOR3314

Regular 4.1




Size	Part #
Ø4.1 x 8mm, Regular Platform One Stage Implant	IOR4108
Ø4.1 x 10mm, Regular Platform One Stage Implant	IOR4110
Ø4.1 x 12mm, Regular Platform One Stage Implant	IOR4112
Ø4.1 x 14mm, Regular Platform One Stage Implant	IOR4114
Ø4.1 x 16mm, Regular Platform One Stage Implant	IOR4116

Regular 4.8




Size	Part #
Ø4.8 x 8mm, Regular Platform One Stage Implant	IOR4808
Ø4.8 x 10mm, Regular Platform One Stage Implant	IOR4810
Ø4.8 x 12mm, Regular Platform One Stage Implant	IOR4812
Ø4.8 x 14mm, Regular Platform One Stage Implant	IOR4814
Ø4.8 x 16mm, Regular Platform One Stage Implant	IOR4816

Wide 4.1



Size	Part #
Ø4.1 x 8mm, Wide Platform One Stage Implant	IOW4108
Ø4.1 x 10mm, Wide Platform One Stage Implant	IOW4110
Ø4.1 x 12mm, Wide Platform One Stage Implant	IOW4112

Wide 4.8




Size	Part #
Ø4.8 x 8mm, Wide Platform One Stage Implant	IOW4808
Ø4.8 x 10mm, Wide Platform One Stage Implant	IOW4810
Ø4.8 x 12mm, Wide Platform One Stage Implant	IOW4812

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Wide 5.6




Size	Part #
Ø5.6 x 8mm, Wide Platform One Stage Implant	ITW5608
Ø5.6 x 10mm, Wide Platform One Stage Implant	ITW5610

Regular 4.8



Size	Part #
Ø4.8 x 6mm, Regular Platform One Stage Implant	IORS

Wide 5.6



Size	Part #
Ø5.6 x 6mm, Wide Platform One Stage Implant	IOWS

Special Cover Screws

Internal Cover Screws



Size	Part #
0mm - Internal Cover Screw, One Stage System	CF0010

1mm Cover Screws



Size	Part #
1mm - Narrow Cover Screw, One Stage System	CN0010
1mm - Regular Cover Screw, One Stage System	CR0010
1mm - Wide Cover Screw, One Stage System	CW0010

3mm Cover Screws



Size	Part #
3mm - Narrow Cover Screw, One Stage System	CN0030
3mm - Regular Cover Screw, One Stage System	CR0030
3mm - Wide Cover Screw, One Stage System	CW0030

Solid Abutments and Restorative Kits

Master Restorative Kits include a solid abutment, one metal analog (for finishing restoration), two resin analogs (for model), snap impression transfer, engaging and non engaging waxing sleeves.

Note: For single crowns use the clear engaging sleeve for undercasting and the white sleeve is used for temporary cap. For splinted restorations use the white sleeve for the undercasting and the clear one for the temporary cap.

All prices are available in the product section of our website. Volume discounts are available for some items.






Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Solid Prosthetics Regular Platform 4.0mm

	Size	Part #
	Master Restorative Kit - Regular 4.0mm	KRT0040
	Abutment, regular 4mm	AR4000
	Snap Impression Transfer, Reg. 4mm	MSR4000
	Abutment Level Analog 4mm	BAR4000
	Waxing/Temporary Sleeves Regular	WIRK4800

Solid Prosthetics Regular Platform 5.5mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KRT0055
	Abutment, regular 5.5mm	AR5500
	Snap Impression Transfer, Reg. 5.5mm	MSR5500
	Abutment Level Analog 5.5mm	BAR5500
	Waxing/Temporary Sleeves Regular	WIRK4800

Solid Prosthetics Regular Platform 7.0mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KRT0070
	Abutment, regular 5.5mm	AR7000
	Snap Impression Transfer, Reg. 5.5mm	MSR7000
	Abutment Level Analog 5.5mm	BAR7000
	Waxing/Temporary Sleeves Regular	WIRK4800

Solid Prosthetics Narrow Platform 4.0mm

	Size	Part #
	Master Restorative Kit - Narrow 4.0mm	KNT0040
	Abutment, regular 5.5mm	AR4000
	Snap Impression Transfer, Narrow 5.5mm	MSN4000
	Abutment Level Analog 4mm	BAN4000
	Waxing/Temporary Sleeves Narrow	WINK4100

Solid Prosthetics Narrow Platform 5.5mm

	Size	Part #
	Master Restorative Kit - Narrow 5.5mm	KNT0055
	Abutment, regular 4mm	AR5500
	Snap Impression Transfer, Reg. 4mm	MSN5500
	Abutment Level Analog 4mm	BAN5500
	Waxing/Temporary Sleeves Narrow	WINK4100

Solid Prosthetics Narrow Platform 7.0mm

	Size	Part #
	Master Restorative Kit - Regular 5.5mm	KNT0070
	Abutment, regular 5.5mm	AR7000
	Snap Impression Transfer, Reg. 5.5mm	MSN7000
	Abutment Level Analog 5.5mm	BASN7000
	Waxing/Temporary Sleeves Regular	WINK4800

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Solid Prosthetics Wide Platform 4.0mm



Size	Part #
Master Restorative Kit - Wide 4.0mm	KWT0040
Abutment, wide 4.0mm	AW4000
Snap Impression Transfer, Wide 4mm	MSW4000
Abutment Level Analog 4.0mm	BAW4000
Waxing/Temporary Sleeves Wide	WIWK6500

Solid Prosthetics Wide Platform 5.5mm



Size	Part #
Master Restorative Kit - Wide 5.5mm	KWT0055
Abutment, wide 5.5mm	AW5500
Snap Impression Transfer, Wide 5.5mm	MSW5500
Abutment Level Analog 5.5mm	BAW5500
Waxing/Temporary Sleeves Wide	WIWK6500

synOcta® Compatible Abutments

Abutments include a fixation screw. Only Compatible with Tissue Level implants. Compatible with Straumann®-ITI® Lifecore® Stage 1® implants. synOcta® compatible. Must be inserted with .048 hex driver. Regular abutments can be used with the regular and narrow platform implants.

synOcta® Compatible Abutments




Size	Part #
Abutment Oct Regular/Narrow	AOR5500
Abutment Oct Wide	AOW5500

Gold Castable Abutments


UCLA abutments are gold alloy. Only Compatible with Tissue Level implants.

Gold UCLA Prosthetics Regular Platform Engaging



Size	Part #
Restorative Kit Gold UCLA - Regular	KAUCLAORG
Gold UCLA Abutment, Regular (including screw)	AUCLAORG
Impression Transfer - Snap Regular	MPOR
Analog - Implant Level Regular	BOR48
Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Regular Platform Engaging



Size	Part #
Restorative Kit Gold UCLA - Regular	KAUCLAORNG
Gold UCLA Abutment Non-Engaging, Regular (including screw)	AUCLORNG
Impression Transfer - Snap Regular	MPOR
Analog - Implant Level Regular	BOR48
Abutment Screw, UCLA, One Stage	STT20UC

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Gold UCLA Prosthetics Wide Platform Engaging

Size	Part #
Restorative Kit Gold UCLA - Wide	KAUCLAOWG
Gold UCLA Abutment, Wide (including screw)	AUCLAOWG
Impression Transfer - Snap Wide	MPOW
Analog – Implant Level Wide	BOW48
Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Wide Platform Non - Engaging

Size	Part #
Restorative Kit Gold UCLA Non - Engaging, Wide	KAUCLAOWNG
Gold UCLA Abutment Non-Engaging, Wide (including screw)	AUCLAOWNG
Impression Transfer - Snap Wide	MPOW
Analog – Implant Level Wide	BOW48
Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Narrow Platform Engaging

Size	Part #
Restorative Kit Gold UCLA - Narrow	KAUCLAONG
Gold UCLA Abutment, Narrow (including screw)	AUCLAONG
Impression Transfer - Snap Narrow	MPON
Analog – Implant Level Narrow	BON48
Abutment Screw, UCLA, One Stage	STT20UC

Gold UCLA Prosthetics Narrow Platform Non - Engaging

Size	Part #
Restorative Kit Gold UCLA - Narrow	KAUCLAONNG
Gold UCLA Abutment, Narrow (including screw)	AUCLONNG
Impression Transfer - Snap Narrow	MPON
Analog – Implant Level Narrow	BON48
Abutment Screw, UCLA, One Stage	STT20UC

Angled Abutments

Only Compatible with Tissue Level implants. Angled Abutments include abutment screw. For assistance in choosing an angled abutment please see the diagram on the web site..

Gold UCLA Prosthetics Narrow Platform Engaging

Size	Part #
Angled 15 degree flat oriented Regular/Narrow	AOA1565A
Angled 15 degree corner oriented Regular/Narrow	AOA1565B



Overdenture Abutments

Master Overdenture Restorative Kit – Chairside includes: overdenture abutment, metal housing, retention inserts, lab processing insert and block out ring.

Master Overdenture Restorative Kit – Laboratory includes: overdenture abutment, metal housing, retention inserts, lab processing insert, block out ring, analog and impression transfer.

For greater flexibility Blue Sky Bio overdenture abutments do not have an internal undercut. Instruments list above are needed for implant and abutment placement. Overdenture abutment can be placed with a .048 hex driver or a Straumann SCS driver.

All prices are available in the product section of our website. Volume discounts are available for some items.




Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann







Super Snap Overdenture Restorative Kits - Chairside Requires No Metal Housing







Super Snap Overdenture Restorative Kits – Laboratory Includes Metal Housing

	Size	Part #
	Super Insert Restorative Kits - Chairside, 1mm	SIRK10
	Super Insert Restorative Kit - Chairside, 2mm	SIRK20
	Super Insert Restorative Kit - Chairside, 3mm	SIRK30

	Size	Part #
	Super Insert Restorative Kit - Laboratory, 1mm	SIRK10-L
	Super Insert Restorative Kit - Laboratory, 2mm	SIRK20-L
	Super Insert Restorative Kit - Laboratory, 3mm	SIRK30-L

Overdenture Products

	Size	Part #
	Complete Flat Snap Restorative Kit	NOPK
	Analog	FSA
	Replacement Inserts Strong Retention (4 Pack)	FSCS
	Replacement Inserts Medium / Light Retention (4 Pack)	FSCG
	Resin Impression transfer (2 Pack) and FSCB (4 Pack)	FSCR
	Processing Inserts No Retention (4 Pack)	FSCB

	Size	Part #
	Super Insert – Strong retention, requires no housing (includes 2 inserts and 2 block out rings) Extreme long lasting retention on parallel implants.	FOPS
	Super Insert – Medium retention, requires no housing (includes 2 inserts and 2 block out rings) Extreme long lasting retention on parallel implants.	FOPG
	Super Insert – Light retention, requires no housing (includes 2 inserts and 2 block out rings) Extreme long lasting retention on parallel implants.	FOPR
	Resin Impression transfer (2 Pack)	FSI
	Block Out Ring (10 Pack)	RSD01
	Insert Placement/Removal Instrument	TFSU

Intraorally Preparable Abutments

Only Compatible with Tissue Level implants. Compatible with Straumann®-ITI® Lifecore® Stage 1® implants. Blue Sky Bio 1-piece angle correction abutments have increased fracture resistance and require intraoral preparation with a carbide bur and copious irrigation. Blue Sky Bio angle correction abutments are not a 2-piece design like ITI-Straumann® angled abutments.

Angle Correction Abutments will fit narrow and regular platform implants.

Must be inserted with a .048 screwdriver head.

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

15° Angle Correction

Size	Part #
6mm - Preparable Abutment Regular/Narrow	AU5541



25° Angle Correction

Size	Part #
6mm - Preparable Abutment Regular/Narrow	AU5548






Waxing Sleeves / Temporary Transfers

Compatible with Straumann®-ITI® Lifecore® Stage 1® implants.

Each package contains one engaging and one non-engaging sleeve, to be used for wax-ups or for temporary transfers..

Waxing Sleeves / Temporary Transfers

Size	Part #
 <p>Clear / Engaging White / Non Engaging</p>	WINK4100
 <p>Clear / Engaging White / Non Engaging</p>	WIRK4800
 <p>Clear / Engaging White / Non Engaging</p>	WIWK6500

Abutment level analogs - Metal

Abutment level analogs - Metal

Size	Part #
4mm Narrow Abutment Analog (BIO One Stage Compatible Only)	BAN4000



5.5mm Narrow Abutment Analog (BIO One Stage Compatible Only)	BAN5500
--	---------



Size	Part #
5.5mm Regular Abutment Analog	BAR5500



7mm Regular Abutment Analog	BAR7000
-----------------------------	---------



Size	Part #
7mm Narrow Abutment Analog (BIO One Stage Compatible Only)	BAN7000



4mm Regular Abutment Analog	BAR4000
-----------------------------	---------



Size	Part #
4mm Wide Abutment Analog	BAW4000



5.5mm Wide Abutment Analog	BAW5500
----------------------------	---------



All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Impression Transfers and Implant Level Analogs

synOcta® compatible.

Must be inserted with .048 hex drive.

Resin Impression Transfers - Snap




Size	Part #
Impression Transfer, Narrow (BIO One Stage Compatible Only)	MPON
Impression Transfer, Regular	MPOR
Impression Transfer, Wide	MPOW

Implant level analogs



Size	Part #
Analog, Oct Narrow (BIO One Stage Compatible Only)	BON41
Analog, Oct Regular	BOR48
Analog, Oct Wide	BOW65

Aluminum Impression Transfers - Closed Tray



Size	Part #
Impression Transfer, Oct Narrow (BIO One Stage Compatible Only)	MON5500
Impression Transfer, Oct Regular	MOR5500
Impression Transfer, Oct Wide	MOW5500

Aluminum Impression Transfers - Open Tray



Size	Part #
Impression Transfer, Oct Narrow (BIO One Stage Compatible Only)	MON5500-OT
Impression Transfer, Oct Regular	MOR5500-OT
Impression Transfer, Oct Wide	MOW5500-OT

Surgical Instruments

Surgical Instruments

















Size	Part #
Torque/Ratchet Wrench	TR-900-01
Thumbwheel/ITI Wrench Adapter	ITI-Adapter
.048 Hex Driver short	SDS048
.048 Hex Driver long	SDL048
J-Implant Insertion Driver short	DT0014

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Surgical Instruments






	Size	Part #
	J-Implant Insertion Driver long	DT0018
	Handpiece Implant Driver short	DMS
	Handpiece Implant Driver long	DML
	Drill Starter, short	DRILLTSS
	Drill Starter, long	DRILLTSL
	Drill 2.2mm, short	DRILLT22S
	Drill 2.2mm, long	DRILLT22L
	Drill 2.8mm, short	DRILLT28S
	Drill 2.8mm, long	DRILLT28L
	Drill 3.5mm, short	DRILLT35S
	Drill 3.5mm, long	DRILLT35L
	Drill 4.3mm, short	DRILLT43S
	Drill 4.3mm, long	DRILLT43L
	Drill 5.0mm, short	DRILLT50S

Safety Stop Surgical Kit



Part #	Part #
Safety Stop Surgical Kit	SSSQ

Safety Stop Kit Individual Components And Accessories

	Size	Part #
	Drill Starter	DRILNSTPS
	Drill Stop 2.2 x 6mm	DRILTSTP2206
	Drill Stop 2.2 x 8mm	DRILTSTP2208
	Drill Stop 2.2 x 10mm	DRILTSTP2210
	Drill Stop 2.2 x 12mm	DRILTSTP2212

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Safety Stop Kit Individual Components And Accessories

	Size	Part #
	Drill Stop 2.8 x 6mm	DRILTSTP2806
	Drill Stop 2.8 x 8mm	DRILTSTP2808
	Drill Stop 2.8 x 10mm	DRILTSTP2810
	Drill Stop 2.8 x 12mm	DRILTSTP2812
	Drill Stop 3.5 x 6mm	DRILTSTP3506
	Drill Stop 3.5 x 8mm	DRILTSTP3508
	Drill Stop 3.5 x 10mm	DRILTSTP3510
	Drill Stop 3.5 x 12mm	DRILTSTP3512
	Drill Stop 4.3 x 6mm	DRILTSTP4306
	Drill Stop 4.3 x 8mm	DRILTSTP4308
	Drill Stop 4.3 x 10mm	DRILTSTP4310
	Drill Stop 4.3 x 12mm	DRILTSTP4312
	Drill Stop Drill Extender	TDE

One - Stage Drill Taps







	Size	Part #
	Tap 3.3mm	TAPT33
	Tap 4.1mm	TAPT41
	Tap 4.8mm	TAPT48
	Tap 5.6mm	TAPT56

One - Stage Guide Pins



	Size	Part #
	Guide Pin 2.2mm	GP0022
	Guide Pin 2.8mm	GP0028
	Guide Pin 3.5mm	GP0035
	Guide Pin 4.3mm	GP0043

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann


Other Surgical Instrumentation











	Size	Part #
	Lindemann Bur	LINDR-S21
	Lance Drillt	LD150
	Trepine Bur	TPHB-B6
	Final Drill for 7.0mm Implant, short	DRILL6653S
	Angled Screw Driver (20Ncm max torque)	SDUGA
	Screw Removal Tool	AFSR

Hammerless Crestal Sinus Lift Kit

Sinus Lift Kit		
	Size	Part #
	Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit
	Addition Guide Pin	GPINMAXNP

BIO | BlueTouch

	Size	Part #
	Touch-Screen Implant Motor (Adjustable Speed 300 – 40,000 RPM; up to 50 Ncm of torque for placement). Does not include handpiece.	BIO-BT

	Size	Part #		Size	Part #
	Saeshin Strong 20:1 Handpiece External Irrigation (button)	ACL(B)-411		Saeshin Strong 1:1 Straight Handpiece External Irrigation	AT-11
	Saeshin Traus 20:1 Non-Optic Handpiece External Irrigation (button)	CRB26XX		Nozzle Clip 1:1 Straight Handpiece	ST-Nozzle
	Anthogyr Impulsion 20:1 Contra Angle (brush type); 50 Ncm tolerance	AHP-85P-I		Autoclavable Tubing Set for BIO BlueTouch	AE-24
	Anthogyr Mont Blanc 20:1 Contra Angle (high efficiency); 80 Ncm tolerance	AHP-85MB-X		Autoclavable Pump Section for AE-24 Tubing Set	AE-24-PUMP
	Anthogyr Mont Blanc 20:1 Contra Angle with Depth Control (high efficiency); 80 Ncm tolerance	AHP-85MB-CX		Synthetic Premium Handpiece Lubricant	LU1011HS
				Lubricant Tool Attachment (necessary for lubricant)	LT1013HS






All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann


Sutures


All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.

Size	Part #	Size	Part #
 PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2	 Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3
 PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2	 Silk 3-0, FS-2 needle, 12/box	S30FS2
 Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3		


Master Bone Expander Kit

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.

Master Bone Expander Kit		Bone Expander Kit	
Size	Part #	Size	Part #
 Master Bone Expander Kit	MBEK	 Bone Expander Kit	BEK

Lerner Ridge Splitting Chisels	
Size	Part #
 3mm Chisel with 2mm blade thickness	RSC32
3mm Chisel with .5mm blade thickness	RSC35
6mm Chisel with 2mm blade thickness	RSC62
6mm Chisel with .5mm blade thickness	RSC65

Guided Sinus Lift Kit

Guided Sinus Lift Kit	
Size	Part #
 Guided Sinus Lift Kit	Sinus Kit

Safety Lateral Sinus Kit

Safety Lateral Sinus Kit	
Size	Part #
 Safety Lateral Sinus Kit	SLSK

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG. Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Guide Tubes

Guide Tubes

Size	Part #	Size	Part #
1.2 Guide Tube (5 Pack) ID=1.2mm OD=2.2mm Height=4mm Lip=1mm of total height	TUBE22-12-4L	1.9 Guide Tube (5 Pack) ID=1.9mm OD=2.9mm Height=4mm Lip=1mm of total height	TUBE29-19-4L
1.4 Guide Tube (5 Pack) ID=1.4mm OD=2.4mm Height=4mm Lip=1mm of total height	TUBE24-14-4L	2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=4mm Lip=1mm of total height	TUBE3-2-4L
1.55 Guide Tube (5 Pack) ID=1.55mm OD=2.55mm Height=7mm	TUBE25-15-4L	2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=5mm Lip=1mm of total height	Tube20
1.6 Guide Tube (5 Pack) ID=1.6mm OD=2.6mm Height=4mm Lip=1mm of total height	TUBE26-16-4L	2.1 Guide Tube (5 Pack) ID=2.1mm OD=3.1mm Height=4mm Lip=1mm of total height	TUBE31-21-4L
1.7 Guide Tube (5 Pack) ID=1.7mm OD=2.7mm Height=4mm Lip=1mm of total height	TUBE27-17-4L	2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-21-5L
2.1 Guide Tube ID=2.1mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-21-5L	2.4 Guide Tube (5 Pack) ID=2.4mm OD=3.4mm Height=4mm Lip=1mm of total height	TUBE34-24-4L
2.1 Guide Tube ID=2.1mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-21-5L	2.4 Guide Tube ID=2.4mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-24-5L
2.11 Guide Tube (5 Pack) ID=2.11mm OD=3.11mm Height=5mm Lip=1mm of total height	Tube211	2.4 Guide Tube ID=2.4mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-24-5L
2.2 Guide Tube (5 Pack) ID=2.2mm OD=3.2mm Height=5mm Lip=1mm of total height	Tube22	2.4 Guide Tube ID=2.4mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-24-5L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Guide Tubes

Size	Part #	Size	Part #
2.3 Guide Tube ID=2.3mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-23-5L	2.6 Guide Tube ID=2.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-26-5L
2.3 Guide Tube ID=2.3mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-23-5L	2.72 Guide Tube ID=2.72mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-272
2.3 Guide Tube ID=2.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-23-5L	2.8 Guide Tube (5 Pack) ID=2.8mm OD=3.8mm Height=5mm Lip=1mm of total height	TUBE280
2.31 Guide Tube (5 Pack) ID=2.31mm OD=3.31mm Height=5mm Lip=1mm of total height	Tube231	2.8 Guide Tube ID=2.8mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-28-5L
2.8 Guide Tube ID=2.8mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-28-5L	3.3 Guide Tube ID=3.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-33-5L
2.9 Guide Tube ID=2.9mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-29-5L	3.5 Guide Tube ID=3.5mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.5 Guide Tube ID=3.5mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.6 Guide Tube ID=3.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-36-5L
2.92 Guide Tube (5 Pack) ID=2.92mm OD=3.92mm Height=5mm Lip=1mm of total height	TUBE292	3.68 Guide Tube (5 Pack) 3.5mm Short/Long, Bio Trilobe ID=3.68mm OD=4.68mm Height=5mm Lip=1mm of total height	TUBE368
3.0 Guide Tube ID=3.0mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-3-5L	3.75 Guide Tube ID=3.75mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-375

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Guide Tubes

Size	Part #	Size	Part #
3.15 Guide Tube ID=3.15mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-315	3.79 Guide Tube ID=3.79mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-379
3.3 Guide Tube ID=3.3mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-33-5L	3.9 Guide Tube ID=3.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-39-5L
4.1 Guide Tube ID=4.1mm OD=5.0mm Height=4mm Lip=1mm of total height	TUBE5-41-4L	4.45 Guide Tube (5 Pack) 4.3mm Short / Long, Bio Trilobe ID=4.45mm OD=5.45mm	TUBE445
4.1 Guide Tube (5 Pack) ID=4.1mm OD=5.1mm Height=3.5mm Lip=1mm of total height	TUBE51-41-35L	4.52 Guide Tube ID=4.52mm OD=5.57mm Height=4mm Lip=1mm of total height	TUBE557-452-4L
4.15 Guide Tube (5 Pack) ID=4.15mm OD=5.15mm Height=5mm Lip=1mm of total height	TUBE415	4.68 Guide Tube ID=4.68mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-468
4.2 Guide Tube ID=4.2mm OD=5.0mm Height=4.5mm Lip=1mm of total height	TUBE5-42-45L	4.7 Guide Tube (5 Pack) ID=4.7mm OD=5.7mm Height=4mm Lip=1mm of total height	TUBE57-47-4L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=4mm Lip=1mm of total height	TUBE525-425-4L	5.0 Guide Tube (5 Pack) ID=5.0mm OD=6.0mm Height=3.5mm Lip=1mm of total height	TUBE6-5-35L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=5mm Lip=1mm of total height	TUBE425	5.0 Guide Tube ID=5.0mm OD=6.2mm Height=4mm Lip=1mm of total height	TUBE62-5-4L
4.25 Guide Tube ID=4.25mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-425	5.02 Guide Tube (5 Pack) ID=5.02mm OD=6.02mm Height=5mm Lip=1mm of total height	TUBE502
4.35 Guide Tube ID=4.35mm OD=5.27mm Height=6mm Lip=1mm of total height	TUBE527-435-6L	5.06 Guide Tube (5 Pack) ID=5.06mm OD=6.02mm Height=3.5mm Lip=1mm of total height	TUBE6-506-35L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Guide Tubes

Size	Part #	Size	Part #
5.08 Guide Tube (5 Pack) ID=5.08mm OD=6.08mm Height=5mm Lip=1mm of total height	TUBE508	5.5 Guide Tube ID=5.5mm OD=6.5mm Height=4mm Lip=1mm of total height	TUBE65-55-4L
5.1 Guide Tube ID=5.1mm OD=6.1mm Height=6mm Lip=1mm of total height	TUBE61-51-6L	5.7 Guide Tube (5 Pack) ID=5.7mm OD=6.7mm Height=4mm Lip=1mm of total height	TUBE67-57-4L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
5.16 Guide Tube (5 Pack) 5.0mm Short, Bio Trilobe ID=5.16mm OD=6.16mm Height=5mm Lip=1mm of total height	TUBE516	6.2 Guide Tube ID=6.2mm OD=6.9mm Height=4mm Lip=1mm of total height	TUBE69-62-4L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=4mm Lip=1mm of total height	TUBE625-525-4L	6.2 Guide Tube (5 Pack) ID=6.2mm OD=7.2mm Height=3.5mm Lip=1mm of total height	TUBE72-62-35L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=5mm Lip=1mm of total height	TUBE525	6.26 Guide Tube (5 Pack) ID=6.26mm OD=7.21mm Height=3.5mm Lip=1mm of total height	TUBE72-626-35L
5.31 Guide Tube ID=5.31mm OD=6.2mm Height=4.5mm Lip=1mm of total height	TUBE62-531-45L	6.3 Guide Tube ID=6.3mm OD=7.3mm Height=6mm Lip=1mm of total height	TUBE73-63-6L
5.35 Guide Tube (5 Pack) ID=5.35mm OD=6.35mm Height=4mm Lip=1mm of total height	TUBE635-535-4L	6.72 Guide Tube ID=6.72mm OD=7.17mm Height=5mm Lip=1mm of total height	TUBE717-672

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Straumann® and Lifecore Biomedical. SinOcta® is a registered trademark of Strauman LG.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Straumann

Terms & Conditions – Important User Information

Blue Sky Bio, LLC - TERMS AND CONDITIONS – Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

NO RETURNS

Blue Sky Bio does not accept returns. Contact us before attempting to return products. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate. All returned products are subject to inspection prior to any credit/refund authorization. Products returned, that do not meet the exceptions listed above, will incur a restocking fee and may be subject to no credit/refund authorization.

A Return Authorization Number (RA#) must be obtained before returning products.

ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 11am CST will ship the same day. Orders placed after 11am CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays, as this will affect shipping delivery dates.

INTERNATIONAL BROKER / TAX TERMS

Blue Sky Bio utilizes Fedex Trade Networks (FTN) to ship and broker all orders from the United States to your door. As broker, FTN pays all taxes, tariffs, and VAT due on shipments. FTN will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FTN of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am – 5pm CST at, **888-446-6724**, or email us anytime at, info@blueskybio.com.

Important Product Notices For Blue Sky Bio Customers

Blue Sky Bio internal taper/oct dental implants are compatible with Straumann® implants, synOcta® system solid abutments and instrumentation. Blue Sky Bio implants have a patent-pending thread design and should not be inserted with mounts from other companies.

Straumann® users require Blue Sky Bio J-implant insertion driver or handpiece driver, and the ITI Wrench Adapter. The J-implant insertion driver eliminates the need for an implant mount. All other instruments are part of the Straumann® Surgical Kit.

Our implants and solid abutment system parts are compatible with Lifecore® Stage-1® implants and components.

With the exceptions listed above no other companies products should be used with our system.

Federal Law restricts these devices to sale by or on order of a licensed dentist or physician.