

Your Partner In Digital Dentistry

BIO | Max



2020

Product Catalog BIO | Max

Nobel Active® Compatible

January 2020

Table of contents

Implants	2
Short Implants	2,3
3.0mm Implants	3
Healing Abutments	3
MPA Custom Abutments	4
KISS Abutments	4
UCLA Gold and Plastic Castable Abutments	4,5
Angled Abutments	5
Multi Unit Abutments	5,6
Standard Temporary Abutment	6
Digital Temporary Abutment	6
Digital Abutment	6
Intraoral Scanning Parts	7
Overdentures - Straight and Angled	7,8
Analogs	8
Impression Transfers	8
Surgical Instruments	8,9

Table of contents

Surgical Instruments	9,10
Safety Stop Surgical Kit	10,11
Fully Guided Keyless Kit	11,12
Hammerless Crestal Sinus Lift Kit	13
BIO BlueTouch	13
Suture	13
Master Bone Expander Kit	14
Guided Sinus Lift Kit	14
Safety Lateral Sinus Kit	14
Guide Tubes	14-18
Terms and Conditions	19

Implants

Packs include BIO | Max implant and 1mm cover screw.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

3.5mm Ø NP platform switched			
	Size	Part #	
	8mm	IJH3508	
	10mm	IJH3510	
	11.5mm	IJH3511	
	13mm	IJH3513	
	16mm	IJH3516	







7.0mm Ø NP platform switched			
	Size	Part #	
	8mm	IJH7008	
	10mm	IJH7010	
	11.5mm	IJH7011	
1			

Short Implants


Kits include BIO | Max implant and 1mm cover screw.


BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

4.3mm Ø NP platform switched			
	Size	Part #	
	6mm	IJH4306	









3.0mm Implants

Kits include BIO | Max implant and 1mm cover screw.

Requires DRILLN26S or DRILLN26L.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.



Healing Abutments

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.




MPA Custom Abutments

MPA Multi Purpose Abutments: Abutments can be used as an implant level impression transfer and the final abutment. Kits include MPA Abutment, Analog and Fixation Screw.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.



KISS Abutments

Fixation screw included

KISS Complete Kits include KISS abutment, screw, metal analog. two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

KISS Impression Kits include two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.





UCLA Gold Castable Abutments

Fixation screw included

Please see the KISS abutment section for whiteMetal abutment screw retained restorations.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.





Plastic Castable Abutments



Abutment NP platform	
Size	Part #
Plastic Non-Engaging UCLA Abutment (including screw	PAUCLJHNPN
Additional Abutment Screw, UCLA, BIO Max	SRJH

Abutment NP platform	
Size	Part #
MPA Plastic Kit Castable Non - Engaging (including abutment, analog and screw)	KPMPAJHNE41
MPA Plastic Abutment Castable Non - Engaging 4.1mm	PMPAJHNE41

Angled Abutments

Fixation screw included.


BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments**.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

Angled Abutments				
		Size	Part #	
		Angled Abutment 15 degree	AAJH15	

Angled Abutments				
		Size	Part #	
		Angled Abutment 30 degree	AAJH30	

Multi Unit Abutments



	Componen	its		
Size	Part #		Size	Part #
Comfort Cap for all Multi Unit Abutments	СМИ		MMUTC Casting Sleeve	PWAXMUTC
Combination Impression Transfer, Temporary Abutment and Castable Sleeve, Plastic - Requires screw SMUS or SMUL	MMUP		Analog, Multi Unit Angled Multi Unit Abutment Delivery Aid	BMU
Combination Impression Transfer/ Temporary Abutment - Titanium Alloy - Requires screw SMUS or SMUL	MMUT		Screw Multi Unit Impression Transfer, long	SMUL
Multiunit Coping, Temporary Abutment Titanium Alloy - Requires screw SMUS or SMUL	MMUTC		Screw, Multi Unit, Castable Abutments, short	SMUS



Standard Temporary Abutments



Size	Part #
Engaging Pickup Abutment for Permanent and Temporary Screwmentable Restorations (includes screw)	ACALJHNP



Size	Part #
Non-Engaging Pickup Abutment for Permanent and Temporary Screwmentable Restorations, (includes screw)	ACALJHNPN

Digital Temporary Abutments

Digital Temporary Abutments						
	Size	Part #		Size	Part #	
	Digital Abutment for temporary abutment, Narrow Platform (includes fixation screw)	ABJHNP		Titanium Digital Abutment - without hex, Narrow Platform (includes abutment screw)	ABJHNPNC	
	Digital Abutment for temporary abutment 1.8mm collar, Narrow Platform (includes fixation screw)	ABJHNP18		Digital Abutment for temporary abutment - non engaging, Narrow Platform Platform (includes abutment screw)	ABJHNPN	
	Digital Abutment for temporary abutment, Mini, Narrow Platform (includes fixation screw)	ABJHNPMI		Angled Screw Driver (20Ncm max torque)	SDUGA	

Digital Abutments

			Digital Abutments		
		D . "	Digital Abadinents	6.	Dout #
	Size	Part #		Size	Part #
James .	Digital Abutment for temporary abutment, Narrow Platform (includes fixation screw)	ABJHNP		Titanium Digital Abutment - without hex, Narrow Platform (includes abutment screw)	ABJHNPNC
	Digital Abutment for temporary abutment 1.8mm collar, Narrow Platform (includes fixation screw)	ABJHNP18		Digital Abutment for temporary abutment - non engaging, Narrow Platform Platform (includes abutment screw)	ABJHNPN
	Digital Abutment for temporary abutment, Mini, Narrow Platform (includes fixation screw)	ABJHNPMI	<u> </u>	Angled Screw Driver (20Ncm max torque)	SDUGA

fixation screw)	BJHNPMI		_		
	Angled	Digital A	butments		
Size		Part #		Size	Part #
Titanium Angled Base Abutment for temporary abutment,Narrow Platfo (includes abutment screw		BJHNP		Titanium Angled Base Abutment for temporary abutment - without hex, Narrow Platform	AABJHNPNC
Titanium Angled Base Abutment for temporary abutment - non engaging, Narrow Platform (includes abutment screw	AAE	BJHNPN		(includes abutment screw)	

Intraoral Scanning Parts

Master Overdenture Restorative Kits Regular - Chairside

	Size	Part #
= (Rosetta Scanpost Extender for Regular Digital Abutments	OPTIBASER
	Regular Ball ScanSnap (required for Exocad or 3Shape)	BALL6
	Rosetta Scanpost Extender for Mini Digital Abutments	OPTIBASEN

Master Overdenture Restorative Kits Regular - Laboratory

5120	
Mini Ball ScanSnap (required for Exocad / 3Shape)	BALL6MI
ScanFlag Scanpost for Regular Digital Abutments	OPTICSTB
ScanFlag Scanpost for Mini	OPTICSTBMI

Overdentures

Master Overdenture Restorative Kits Regular - Chairside

Size	Part #
Master Overdenture Restorative Kit - Chairside, regular x 3.0mm height	MORKJH10
Master Overdenture Restorative Kit - Chairside, regular x 3.0mm height	MORKJH30
Master Overdenture Restorative Kit - Chairside, regular x 5.0mm height	MORKJH50

Master Overdenture Restorative Kits Regular - Laboratory

Master Overdenture Restorative Rits Regular - Laboratory				
	Size	Part #		
	Master Overdenture Restorative Kit - Laboratory, regular x 3.0mm height	MORKJH10-L		
	Master Overdenture Restorative Kit - Laboratory, regular x 3.0mm height	MORKJH30-L		
	Master Overdenture Restorative Kit - Laboratory, regular x 5.0mm height	MORKJH30-L		

Angled Overdentures

Angled Overdenture Restorative Kits Regular - Chairside

Size	Part #
Angled Overdenture Restorative Kit , 17 degrees - Chairside	ARKJH17

Angled Overdenture Restorative Kits Regular - Laboratory

Size	Part #
Angled Overdenture Restorative Kit, 17 degrees - Laboratory	ARKJH17-L

Angled Super Snap Overdenture Restorative Kits - Chairside Requires No Metal Housing

Kit, 17 degrees - Chairside

Size	Part #
Angled Super Overdenture Restorative Kit , 17 degrees - Chairside	ASRKJH17
Angled Super Overdenture Restorative	

Angled Super Snap Overdenture Restorative Kits - Laboratory Requires No Metal Housing

Angled Super Overdenture Restorative Kit , 17 degrees - Laboratory	ASRKJH17-L
Angled Super Overdenture Restorative Kit , 17 degrees - Laboratory	ASRKJH30-L

Part #

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare.

ASRKJH30

Overdenture Products Size Part # Size Part # Super Insert - Light retention, requires no housing(includes 2 inserts and 2 **FOPR** Complete Flat Snap Restorative Kit NOPK block out rings) Extreme long lasting retention on parallel implants. FSI Resin Impression transfer (2 Pack) **FSA** Analog Replacement Inserts Strong **FSCS** Retention (4 Pack) **HGFS** Titanium Housing Replacement Inserts Medium / Light **FSCG** Retention (4 Pack) RSD01 Resin Impression transfer (2 Pack) Block Out Ring (10 Pack) **FSCR** and FSCB (4 Pack) Processing Inserts No Retention **FSCB TFSU** Insert Placement/Removal Instrument (4 Pack) Super Insert – Strong retention, requires no housing (includes 2 inserts and 2 **FOPS** block out rings) Extreme long lasting **TFSU** Insert Placement/Removal Instrument retention on parallel implants. Super Insert - Medium retention, requires **FOPG** no housing(includes 2 inserts and 2

Analogs



Impression Transfers

block out rings)Extreme long lasting retention on parallel implants.



Surgical Instruments






Other Surgical Instrumentation




Other Surgical Instrumentation			
Size	Part #		
Final Drill for 7.0mm Implant, short	DRILL6653S		
Final Drill for 7.0mm Implant, short	DRILL6653S		
Angled Screw Driver (20Ncm max torque)	SDUGA		
Screw Removal Tool	AFSR		

Safety Stop Surgical Kit



	Safety Stop Kit Individual Components And Accessories		
	Size	Part #	
Ø3.9 -	Drill Stop 3.9 x 8mm	DRILNSTP3908	
Ø3.9°=	Drill Stop 3.9 x 10mm	DRILNSTP3910	
Ø3.9 ₋	Drill Stop 3.9 x 11.5mm	DRILNSTP3911	
Ø3.9° -	Drill Stop 3.9 x 13mm	DRILNSTP3913	
Ø4.6 -	Drill Stop 4.6 x 6mm	DRILNSTP4606	
Ø4.6 -	Drill Stop 4.6 x 8mm	DRILNSTP4608	
Ø4.6 -J	Drill Stop 4.6 x 10mm	DRILNSTP4610	
Ø4.6 =	Drill Stop 4.6 x 11.5mm	DRILNSTP4611	
Ø4.6 -	Drill Stop 4.6 x 13mm	DRILNSTP4613	
	Drill Stop Drill Extender	TDE	

Fully Guided Keyless Kit



All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

Guided Kit Individual Components And Accessories Size Part # Guide Drill 2.5 x 8mm DRILG2508 Guide Drill 2.5 x 10mm DRILG2510 Guide Drill 2.5 x 11.5mm DRILG2011 Guide Drill 2.5 x 13mm **DRILG2513** 23.0 Guide Drill 3.0 x 6mm DRILG3006 Guide Drill 3.0 x 8mm DRILG3008 Ø3.0 = Guide Drill 3.0 x 10mm **DRILG3010 03.0** = Guide Drill 3.0 x 11.5mm **DRILG3011** Ø3.0 = Guide Drill 3.0 x 13mm DRILG3013 Guide Drill 3.9 x 6mm Ø3.9 = DRILG3906 Ø3.9 = Guide Drill 3.9 x 8mm **DRILG3908** 03.9 = Guide Drill 3.9 x 10mm **DRILG3910** Guide Drill 3.9 x 11.5mm 03.9 = **DRILG3911** Guide Drill 3.9 x 13mm **DRILG3913** Guide Drill 4.6 x 6mm Ø4.6 = DRILG4606 Guide Drill 4.6 x 8mm Ø4.6 **DRILG4608** Guide Drill 4.6 x 10mm **DRILG4610** 04.6 Guide Drill 4.6 x 11.5mm DRILG4611 Guide Drill 4.6 x 13mm 04.6 **DRILG4613** Bone Crest Flattening Drill **DRILGFLAT Guided Implant Driver** DTGJHNP85 8.5mm Offset Guided Implant Driver **DTGJHNPVR** Variable Offset PIN Guide Fixation Pin (5 Pack)

Hammerless Crestal Sinus Lift Kit

	Sinus Lift Kit	
	Size	Part #
Vinou	Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit
	Addition Guide Pin	GPINMAXNP

BIO | BlueTouch



Size	Part #
Touch-Screen Implant Motor (Adjustable Speed 300 – 40,000 RPM; up to 50 Ncm of torque for placement). Does not include handpiece.	BIO-BT





Sutures

All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.

	Size	Part #
PILITURE IN STREET	PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2
Description of the second of t	PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2
Sep 5 - 0 P Sep 5 - 0	Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3

	Size	Part #
Lists 5-0 EP 1 1900' 10c EP 1 1900' 10c A 20 2000' 10c A 30 2000' 10c A 3	Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3
SADE	Silk 3-0, FS-2 needle, 12/box	S30FS2

Master Bone Expander Kit

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.







Guided Sinus Lift Kit



Safety Lateral Sinus Kit

	Safety Lateral Sinus Kit	
Soliets Lineared Stones Alt	Size	Part #
	Safety Lateral Sinus Kit	SLSK

Guide Tubes

		Guide Tubes		
Size	Part #		Size	Part #
1.2 Guide Tube (5 Pack) ID=1.2mm OD=2.2mm Height=4mm Lip=1mm of total height	TUBE22-12-4L		1.9 Guide Tube (5 Pack) ID=1.9mm OD=2.9mm Height=4mm Lip=1mm of total height	TUBE29-19-4L
1.4 Guide Tube (5 Pack) ID=1.4mm OD=2.4mm Height=4mm Lip=1mm of total height	TUBE24-14-4L		2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=4mm Lip=1mm of total height	TUBE3-2-4L

Size	Part #	Size	Part #
1.55 Guide Tube (5 Pack) ID=1.55mm OD=2.55mm Height=7mm		2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=5mm Lip=1mm of total height	Tube20
1.6 Guide Tube (5 Pack) ID=1.6mm OD=2.6mm Height=4mm Lip=1mm of total height	TUBE26-16-4L	2.1 Guide Tube (5 Pack) ID=2.1mm OD=3.1mm Height=4mm Lip=1mm of total height	TUBE31-21-4L
1.7 Guide Tube (5 Pack) ID=1.7mm OD=2.7mm Height=4mm Lip=1mm of total height	TUBE27-17-4L	2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-21-5L
2.1 Guide Tube ID=2.1mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-21-5L	2.4 Guide Tube (5 Pack) ID=2.4mm OD=3.4mm Height=4mm Lip=1mm of total height	TUBE34-24-4L
2.1 Guide Tube ID=2.1mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-21-5L	2.4 Guide Tube ID=2.4mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-24-5L
2.11 Guide Tube (5 Pack) ID=2.11mm OD=3.11mm Height=5mm Lip=1mm of total height	Tube211	2.4 Guide Tube ID=2.4mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-24-5L
2.2 Guide Tube (5 Pack) ID=2.2mm OD=3.2mm Height=5mm Lip=1mm of total height	Tube22	2.4 Guide Tube ID=2.4mm OD=5.0mm Height=5mm Lip=1mm of total heigh	TUBE5-24-5L
2.3 Guide Tube ID=2.3mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-23-5L	2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE4-26-5L
2.3 Guide Tube ID=2.3mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-23-5L	2.72 Guide Tube ID=2.72mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-272
2.3 Guide Tube ID=2.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-23-5L	2.8 Guide Tube (5 Pack) ID=2.8mm OD=3.8mm Height=5mm Lip=1mm of total height	TUBE280

Size	Part #	Size	Part #
2.31 Guide Tube (5 Pack) ID=2.31mm OD=3.31mm Height=5mm Lip=1mm of total height	Tube231	2.8 Guide Tube ID=2.8mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-28-5L
2.8 Guide Tube ID=2.8mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-28-5L	3.3 Guide Tube ID=3.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-33-5L
2.9 Guide Tube ID=2.9mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-29-5L	3.5 Guide Tube ID=3.5mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.5 Guide Tube ID=3.5mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.6 Guide Tube ID=3.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-36-5L
2.92 Guide Tube (5 Pack) ID=2.92mm OD=3.92mm Height=5mm Lip=1mm of total height	TUBE292	3.68 Guide Tube (5 Pack) 3.5mm Short/Long, Bio Trilobe ID=3.68mm OD=4.68mm Height=5mm Lip=1mm of total height	TUBE368
3.0 Guide Tube ID=3.0mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-3-5L	3.75 Guide Tube ID=3.75mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-375
3.15 Guide Tube ID=3.15mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-315	3.79 Guide Tube ID=3.79mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-379
3.3 Guide Tube ID=3.3mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-33-5L	3.9 Guide Tube ID=3.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-39-5L
4.1 Guide Tube ID=4.1mm OD=5.0mm Height=4mm Lip=1mm of total height	TUBE5-41-4L	4.45 Guide Tube (5 Pack) 4.3mm Short / Long, Bio Trilobe ID=4.45mm OD=5.45mm	TUBE445

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

** BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

Size	Part #	Size	Part #	
4.1 Guide Tube (5 Pack) ID=4.1mm OD=5.1mm Height=3.5mm Lip=1mm of total height	TUBE51-41-35L	4.52 Guide Tube ID=4.52mm OD=5.57mm Height=4mm Lip=1mm of total height	TUBE557-452-	
4.15 Guide Tube (5 Pack) ID=4.15mm OD=5.15mm Height=5mm Lip=1mm of total height	TUBE415	4.68 Guide Tube ID=4.68mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-468	
4.2 Guide Tube ID=4.2mm OD=5.0mm Height=4.5mm Lip=1mm of total height	TUBE5-42-45L	4.7 Guide Tube (5 Pack) ID=4.7mm OD=5.7mm Height=4mm Lip=1mm of total height	TUBE57-47-4L	
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=4mm Lip=1mm of total height	TUBE525-425-4L	5.0 Guide Tube (5 Pack) ID=5.0mm OD=6.0mm Height=3.5mm Lip=1mm of total height	TUBE6-5-35L	
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=5mm Lip=1mm of total height	TUBE425	5.0 Guide Tube ID=5.0mm OD=6.2mm Height=4mm Lip=1mm of total height	TUBE62-5-4L	
4.25 Guide Tube ID=4.25mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-425	5.02 Guide Tube (5 Pack) ID=5.02mm OD=6.02mm Height=5mm Lip=1mm of total height	TUBE502	
4.35 Guide Tube ID=4.35mm OD=5.27mm Height=6mm Lip=1mm of total height	TUBE527-435-6L	5.06 Guide Tube (5 Pack) ID=5.06mm OD=6.02mm Height=3.5mm Lip=1mm of total height	TUBE6-506-35	
5.08 Guide Tube (5 Pack) ID=5.08mm OD=6.08mm Height=5mm Lip=1mm of total height	TUBE508	5.5 Guide Tube ID=5.5mm OD=6.5mm Height=4mm Lip=1mm of total height	TUBE65-55-4L	
5.1 Guide Tube ID=5.1mm OD=6.1mm Height=6mm Lip=1mm of total height	TUBE61-51-6L	5.7 Guide Tube (5 Pack) ID=5.7mm OD=6.7mm Height=4mm Lip=1mm of total height	TUBE67-57-4L	
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-3.	

Size	Part #	Size	Part #
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-6
5.16 Guide Tube (5 Pack) 5.0mm Short, Bio Trilobe ID=5.16mm OD=6.16mm Height=5mm Lip=1mm of total height	TUBE516	6.2 Guide Tube ID=6.2mm OD=6.9mm Height=4mm Lip=1mm of total height	TUBE69-6
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=4mm Lip=1mm of total height	TUBE625-525-4L	6.2 Guide Tube (5 Pack) ID=6.2mm OD=7.2mm Height=3.5mm Lip=1mm of total height	TUBE72-€
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=5mm Lip=1mm of total height	TUBE525	6.26 Guide Tube (5 Pack) ID=6.26mm OD=7.21mm Height=3.5mm Lip=1mm of total height	TUBE72-6
5.31 Guide Tube ID=5.31mm OD=6.2mm Height=4.5mm Lip=1mm of total height	TUBE62-531-45L	6.3 Guide Tube ID=6.3mm OD=7.3mm Height=6mm Lip=1mm of total height	TUBE73-
5.35 Guide Tube (5 Pack) ID=5.35mm OD=6.35mm Height=4mm Lip=1mm of total height	TUBE635-535-4L	6.72 Guide Tube ID=6.72mm OD=7.17mm Height=5mm Lip=1mm of total height	TUBE71

Terms and Conditions

Blue Sky Bio, LLC - TERMS AND CONDITIONS – Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

RETURNS & EXCHANGES

Blue Sky Bio will accept returns or exchanges within 30 days of original sales receipt. Products returned or exchanged between 31 and 120 days of original sales receipt are subject to a 15% restocking fee. No products can be returned or exchanged after 120 days from date of original sales receipt. Prior to return or exchange please contact us for a Return Authorization Number (RA#). All returned or exchanged products are subject to inspection. Products that have been opened, damaged or marked upon are not eligible for return or exchange. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate.

ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 1:30pm CST will ship the same day. Orders placed after 1:30pm CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays and shipping delays due to weather..

INTERNATIONAL BROKER / TAX TERMS

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am – 5pm CST at, **888-446-6724**, or email us anytime at, **info@blueskybio.com**.

IMPORTANT PRODUCT NOTICES FOR BLUE SKY BIO CUSTOMERS

BIO | Max is NobelActive® compatible.

BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

Blue Sky Bio is not affiliated with Nobel Biocare®. NobelActive® is a registered trademark of Nobel Biocare.