


Enbios
World's fastest autoclave.
Vacuum+drying.



Bone dry pouches in no time.



Enbio can sterilize a load of wrapped instruments in just 15 minutes.

Full drying means that pouches are bone dry. Not somewhat dry.

In case of emergency, when you need an instrument sterilized for immediate use, FAST process does the job in just 7 minutes.

15 min

7 min

4 full cycles with drying per hour.



Enbio S is 6 times smaller and 3 times lighter than a typical vacuum autoclave.
But when it comes to work efficiency, it keeps up the pace.

	Pouches / h	Instruments / h
Enbio ^s	32	50



The new definition of sterilization speed.



Time is a relative value. It flows differently under different circumstances.


But when you need to have your instruments ready, you need them, well, FAST.

While other autoclaves might test your patience, Enbio bends the laws of physics to cut the waiting to the invisible.



Typical autoclave





Enbios. Speedmaster.

Vacuum
Pouches
+drying

FDA Cleared

Enbio S is compact to fit even a confined space. Unbelievably short process times will save your procedure in emergency, or allow sterilization of the instruments as you go, between patients.

Unwrapped FAST / 134°C

7 min

Vacuum / 134°C

Wrapped Vacuum / 121°C

15 min

Wrapped

45 min







110V/60Hz

3,23 kW

Unit weight

33 lb

Capacity

Max. packages amount (circa)

Max. recommended batch weight

Process chamber volume

8

1,1 lb

2,7 l

Pressure & temperature

Working pressure 2,1 Bar
Max. pressure 2,3 Bar
Max. process tempreature 137°C

Dimensions


Width 9.8 in Height 6.3 in Depth 22 in

Electrical

Power supply Installed power Max. power consumption Misc


Ingress protection level Data archiving Water conductivity* IP20 USB stick < 15 µS/cm

Enbio user's manual




Enbio Care

Perfect technology meets perfect care.
All Enbio devices are covered with world's unique warranty and care program, with backup device delivered to you for the time of repair.



20.000 units in 4 years.

In just 4 years, Enbio's worldwide community has grown to 20,000 users, who sterilize instruments in our autoclaves, for the complete safety of their customers and patients.













Enbio Group AGEichengasse 3, 4702 Oensingen Switzerland

FDA Establishment Registration Number: 3016777357

US Distribution Center Enbio Corp. 3535 Gravel Springs Rd, Ext. 205 Buford, GA 30519