


Your Partner In Digital Dentistry

# **BIO** Internal Hex





# **Product Catalog** BIO Internal Hex

Zimmer Tapered Screw Vent® Compatible

January 2020

Table of contents	
Implant Packs	2
Practice Building Kits	2,3
Healing Abutments	3
MPA Custom Abutments	4
Angled Abutments	4
KISS Abutments	4,5
Gold Castable Abutments	5
Plastic Castable Abutments	5
Multi Unit Abutments	6
Temporary Abutments	6
Digital Abutments	7
Intraoral Scanning Parts	7
Overdenture Abutments	7
Angled Overdenture Abutments	8
Analog	9
Impression Transfers	9
Surgical Instruments	9,10

Internal Hex Drill Taps	10
Safety Stop Surgical Kit	11,12
Fully Guided Keyless Kit	12,13
Hammerless Crestal Sinus Lift Kit	14
BIO   BlueTouch	14
Suture	14
Master Bone Expander Kit	15
Guided Sinus Lift Kit	15
Safety Lateral Sinus Kit	15
Guide Tubes	15-18
Terms & Conditions	19

### **Implant Packs**

• Packs include Internal Hex implant, 1mm cover screw and a 3mm healing abutment

Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.

3.7mm Implant - 3.5mm Platform				
		Size	Part #	
	ŢĪ	8mm	IH3708	
		10mm	IH3710	
9 _		11.5mm	IH3711	
		13mm	IH3713	
		16mm	IH3716	

4.7mm Implant - 3.5mm Platform Switched				
	Size	Part #		
	6mm	IHPS4706		
	8mm	IHPS4708		
	10mm	IHPS4710		
	11.5mm	IHPS4711		
	13mm	IHPS4713		
	16mm	IHPS4716		

5.7mm Implant - 4.5mm Platform Switched				
	Size	Part #		
	6mm	IH5708		
	8mm	IH5708		
	10mm	IH5710		
	11.5mm	IH5711		
	13mm	IH5713		
	16mm	IH5716		

4.1mm Implant - 3.5mm Platform				
- E - T	Size	Part #		
	8mm	IH4108		
	10mm	IH4110		
	11.5mm	IH4111		
	13mm	IH4113		
	16mm	IH4116		

4.7mm Implant - 4.5mm Platform				
	Size	Part #		
	8mm	IH4708		
	10mm	IH4710		
	11.5mm	IH4711		
	13mm	IH4713		
$\bigcirc$	16mm	IH4716		



### **Practice Building Kits**

• Kits include Internal Hex implant, 1mm cover screw, 3mm abutment, final abutment/implant level impression coping, and implant level analog

• MPA Multi Purpose Abutments: Abutments can be used as an implant level impression transfer and the final abutment. Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.



4.1mm Practice Building Kit - 3.5mm Platform					
		Size	Part #		
	5	8mm	KZPACTA4108		
		10mm	KZPACTA4110		
	( )	11.5mm	KZPACTA4111		
		13mm	KZPACTA4113		
		16mm	KZPACTA4116		

	Size	Part #			Size	Part #
E	6mm	KZPACTAPS4706		- 8	8mm	KZPACTA4708
11	8mm	KZPACTAPS4708			10mm	KZPACTA4710
	10mm	KZPACTAPS4710			11.5mm	KZPACTA471
	11.5mm	KZPACTAPS4711			13mm	KZPACTA4713
	13mm	KZPACTAPS4713			16mm	KZPACTA4716
ctice I	Building Kit - 4.5m Size	-	5.7mm	Practice Bui		5mm Platform Part #
e	Building Kit - 4.5m Size 6mm	n Platform Part # KZPACTA5706	5.7mm	Practice Bui	ilding Kit - 4. Size 6mm	5mm Platform Part <del>1</del> KZPACTA57
	Size	Part #	5.7mm	Practice Bui	Size	Part <del>1</del>
ce I	Size 6mm	Part # KZPACTA5706	5.7mm	Practice Bui	Size 6mm	Part # KZPACTA57
	Size 6mm 8mm	Part # KZPACTA5706 KZPACTA5708	5.7mm	Practice Bui	Size 6mm 8mm	RZPACTA5
3	Size 6mm 8mm 10mm	Part # KZPACTA5706 KZPACTA5708 KZPACTA5710	5.7mm	Practice Bui	Size 6mm 8mm 10mm	Part # KZPACTA57 KZPACTA57 KZPACTA57

### **Healing Abutments and Cover Screws**

Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent<sup>®</sup> implant system.

3.5 Platform Healing Abutment	
Size	Part #
3.5 Platform - 3mm height	HH3530
3.5 Platform - 5mm height	HH3550
3.5 Platform Healing Abutment	
Size	Part #
7mm 3.5 Platform - 3mm height	HH354730
7mm 3.5 Platform - 5mm height	HH354750
7mm 3.5 Platform - 5mm height	HH356550
7mm 3.5 Platform - 5mm height	HH356550
over Screw 3.5mm Platform	CH3500

### **MPA Custom Abutments**

• Kits include MPA Abutment, Analog, and Fixation Screw.

• MPA Multi Purpose Abutments: Abutments can be used as an implant level impression transfer and the final abutment. Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.





	4.5 Platform Abutment	
	Size	Part #
7 1	4.5 MPA Abutment Analog Kit	KMPAIH4547
	4.5 MPA Abutment	MPAIH4547
	4.5 Implant Level Analog	BIH45

	4	4.5 Platform wide Adutment	
		Size	Part #
1		4.5 MPA Wide Abutment Analog Kit	KMPAIH4565
		4.5 MPA Wide Abutment	MPAIH4565
		4.5 Implant Level Analog	BIH45

# **Angled Abutments**

• Fixation screw included

Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.

3.5 Platform Angled Abutment				
H. H.	Size	Part #		
	Angled Abutment 15 degree	AAH3515		
	Angled Abutment 25 degree	AAH3525		

4.5 Platform Angled Abutment				
Au Au	Size	Part #		
	Angled Abutment 15 degree	AAH4515		
	Angled Abutment 25 degree	AAH4525		
	Angled Abutment Wide 15 degree	AAH456515		

# **KISS Abutments**

- Complete Kits include KISS abutment, screw, metal analog, two resin analogs, impression coping, engaging and non-engaging waxing sleeves and temporary sleeve.
- Impression Kits include two resin analogs, impression coping, engaging and non-engaging waxing sleeves and temporary sleeve.

Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.

	3.5 Platform Regular KISS				4.5 Platform Regular KISS	
	Size	Part #			Size	Part #
Æ	3.5 Regular KISS Abutment 1mm collar ASIH35R01		4.5 Regular KISS Abutment 1mm collar	ASIH45R01		
				4.5 Regular KISS Abutment 2mm collar	ASIH45R02	
					4.5 Regular KISS Abutment 3mm collar	ASIH45R03
		4.5 Pla	tform Wide KI	ISS		
		Size			Part #	
		4.5 Wide	KISS Abutment 2mn	n collar	ASIH45W02	

# **Gold Castable Abutments**

Fixation screw included

• UCLA Abutment contains an implant level impression feature. No separate implant level impression coping is required. Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent<sup>®</sup> implant system.

3	.5 Platform Gold UCLA Engaging		4.5	5 Platform Gold UCLA Engaging	
	Size	Part #		Size	Part #
	3.5 Gold UCLA Abutment	AUCLAIH35G		4.5 Gold UCLA Abutment	AUCLAIH45G
	3.5 UCLA Abutment	BIH35		4.5 UCLA Abutment	BIH45
	3.5 Implant Level Analog	SIH01A		4.5 Implant Level Analog	SIH01A

3.5 Platform Gold Non-Engaging UCLA					
		in	Size	Part #	
			3.5 Gold Non-Engaging UCLA	AUCLAIH35NG	
	μ.		3.5 UCLA Abutment	BIH35	
100			3.5 Implant Level Analog	SIH01A	

4.5 Flatform Gold OCLA Engaging						
			Size	Part #		
			4.5 Gold UCLA Abutment	AUCLAIH45G		
			4.5 UCLA Abutment	BIH45		
			4.5 Implant Level Analog	SIH01A		

#### 4.5 Platform Gold Non-Engaging UCLA

Size	Part #
4.5 Gold Non-Engaging UCLA	AUCLAIH45NG
4.5 UCLA Abutment	BIH45
4.5 Implant Level Analog	SIH01A

### **Plastic Castable Abutments**



### **Multi Unit Abutments**

	3.5 Platform						
	Size	Part #					
Ŷ	Abutment Multi Unit Straight, 3.5 platform, Intrenal Hex, 1.5mm cuff - Requires .080 Hex Driver	AMUSIH3515					
Ŷ	Abutment Multi Unit Straight, 3.5 platform, Internal Hex, 3.0mm cuff - Requires .080 Hex Driver	AMUSIH3530					
Ŷ	Abutment Multi Unit Straight, 3.5 platform, Internal Hex, 5.0mm cuff - Requires .080 Hex Driver	AMUSIH3550					
	Angled Abutment Multi Unit, 3.5 platform, 17 degree, Internal Hex, includes abutment screw	AAMUIH3517					
	Angled Abutment Multi Unit, 3.5 platform, 30 degree, Internal Hex, includes abutment screw	AAMUIH3530					

4.5 Platform	
Size	Part #
Abutment Multi Unit Straight, 4.5 platform, Internal Hex, 1.5mm cuff - Requires .080 Hex Driver	AMUSIH4515
Abutment Multi Unit Straight, 4.5 platform, Internal Hex, 3.0mm cuff - Requires .080 Hex Driver	AMUSIH4530
Abutment Multi Unit Straight, 4.5 platform, Internal Hex, 5.0mm cuff - Requires .080 Hex Driver	AMUSIH4550
Angled Abutment Multi Unit, 4.5 platform, 17 degree, Internal Hex, includes abutment screw	AAMUIH4517
Angled Abutment Multi Unit, 4.5 platform, 30 degree, Internal Hex, includes abutment screw	AAMUIH4530


	Components									
	Size	Part #		Size	Part #					
	Comfort Cap for all Multi Unit Abutments	CMU		MMUTC Casting Sleeve	PWAXMUTC					
	Combination Impression Transfer, Temporary Abutment and Castable Sleeve, Plastic - Requires screw SMUS	MMUP		Analog, Multi Unit	BMU					
	or SMUL Combination Impression Transfer/			Angled Multi Unit Abutment Delivery Aid	DSMU					
B	Temporary Abutment - Titanium Alloy - Requires screw SMUS or SMUL	MMUT		Screw Multi Unit Impression Transfer, long	SMUL					
	Multiunit Coping, Temporary Abutment Titanium Alloy - Requires screw SMUS or SMUL	MMUTC		Screw, Multi Unit, Castable Abutments, short	SMUS					
	Multi Unit Titanium Coping Thimble	MUT		.080 Hex Driver for straight Multi Unit Abutments	HFD					
	MMUTC Spacer Sleeve	PSPACEMUTC								

### **Temporary Abutment**

	3.5 Platform Abutmer	it		4.5 Platform Abutme	nt
11	Size	Part #	1	Size	Part #
	Digital Abutment for temporary abutment (includes abutment screw)	ABIH35	hunu	Digital Abutment for temporary abutment (includes abutment screw)	ABIH45
ŧ	Additional Abutment Screw	SIH01A		Additional Abutment Screw	SIH01A

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Zimmer Dental<sup>®</sup>.

### **Digital Abutment**





#### 4.5 Platform Abutment

Size	Part #
Digital Abutment for abutment and crowns for zirconia and emax (includes fixation screw)	ABIH45
Additional Abutment Screw	SIH01A

### **Intraoral Scanning Parts**

Master Overdenture Restorative Kits Regular - Chairside Master Overdenture Restorative Kits Regular - Laboratory Part # Part # Size Size Rosetta Scanpost Extender for Regular Mini Ball ScanSnap OPTIBASER BALL6MI **Digital Abutments** (required for Exocad / 3Shape) Regular Ball ScanSnap ScanFlag Scanpost for Regular OPTICSTB BALL6 (required for Exocad or 3Shape) **Digital Abutments** Rosetta Scanpost Extender for Mini ScanFlag Scanpost for Mini OPTICSTBMI OPTIBASEN **Digital Abutments Digital Abutments** 

#### **Overdenture Abutments**

Super Snap Overdenture Restorative Kit – Chairside includes: overdenture abutment, retention inserts and block out ring. Super Snap Overdenture Restorative Kit – Laboratory includes: overdenture abutment, retention inserts, block out ring, analog and impression transfer.

Super Snap	<b>Overdenture Restorative Kits</b>	3.5- Chairside	Super Snap O	verdenture Restorative Kits 3.	5 - Laboratory
	Size	Part #		Size	Part #
	Super Insert Restorative Kits - Chairside, 3mm	SIRKIH3530		Super Insert Restorative Kits - Laboratory, 3mm	SIRKIH3530-L
	Super Insert Restorative Kits - Chairside, 5mm	SIRKIH3550		Super Insert Restorative Kits - Laboratory, 5mm	SIRKIH3550-L
Super Snap	Overdenture Restorative Kits	4.5 - Chairside	Super Snap Ov	verdenture Restorative Kits 4.	5 - Laboratory
	Size	Part #		Size	Part #
<b>-</b>	Super Insert Restorative Kits -	SIRKIH4530	<b></b>	Super Insert Restorative Kits -	
◎ 🚆 '	Chairside, 3mm	SII ((III 14000		Laboratory, 3mm	SIRKIH4530-L

# **Angled Overdenture Abutments**

Kit 4.5 Platform, 17 degrees ASRKIH4530

r Snap	Overdenture Restorative Kits 3.	5- Chairside	Super Snap O	verdenture Restorative Kits 3.5 -	Laboratory
	Size	Part #		Size	Part #
	Angled Super Overdenture Restorative Kit 3.5 Platform, 17 degrees	ASRKIH3517	₀⊨∎∮	Angled Super Overdenture Restorative Kit 3.5 Platform, 30 degrees	ASRKIH3517-
	Angled Super Overdenture Restorative Kit 3.5 Platform, 30 degrees	ASRKIH3530		Angled Super Overdenture Restorative Kit 3.5 Platform, 30 degrees	ASRKIH3530-
p (	Overdenture Restorative Kits 4.5	- Chairside	Super Snap Ov	verdenture Restorative Kits 4.5 -	Laboratory
	Size	Part #		Size	Part #
	Angled Super Overdenture Restorative Kit 4.5 Platform, 17 degrees	ASRKIH4517	●H <b></b> 4	Angled Super Overdenture Restorative Kit 4.5 Platform, 17 degrees	ASRKIH4517-
2	Angled Super Overdenture Restorative		•••== d	Angled Super Overdenture Restorative	

**Overdenture Products** 

0

Kit 4.5 Platform, 17 degrees

ASRKIH4530-L

Size	Part #		Size
Complete Flat Snap Restorative Kit	NOPK		Super Insert – Light retention, requires no housing(includes 2 inserts and 2 block out rings) <b>Extreme long lasting</b> <b>retention on parallel implants</b> .
Analog	FSA		Resin Impression transfer (2 Pack)
Replacement Inserts Strong Retention (4 Pack)	FSCS		The sines Henrice
Replacement Inserts Medium / Light Retention (4 Pack)	FSCG		Titanium Housing
Resin Impression transfer (2 Pack) and FSCB (4 Pack)	FSCR	0	Block Out Ring (10 Pack)
Processing Inserts No Retention (4 Pack)	FSCB		Insert Placement/Removal Instrument
Super Insert – Strong retention, requires no housing (includes 2 inserts and 2 block out rings) <b>Extreme long lasting</b> retention on parallel implants.	FOPS		Insert Placement/Removal Instrument
Super Insert – Medium retention, requires no housing(includes 2 inserts and 2 block out rings)Extreme long lasting retention on parallel implants.	FOPG		

### Analogs

Implant Level Analogs				
Size	Part #		Size	Part #
3.5 Implant Level Analog	BIH35		4.5 Implant Level Analog	BiH45

### **Impression Transfers**

Implant Level Impression Transfers					
	Size	Part #		Size	Part #
	3.5 Implant Level Impression Transfer	MIH3500		4.5 Implant Impression Transfer	MIH4500

### **Surgical Instruments**

Blue Sky Bio Internal Hex implants are fully compatible with Zimmer Dental's Tapered Screw Vent® implant system.

	Surgical Instruments			
	Size	Part #		
D P BhueSkyBio	Torque/Ratchet Wrench	TR-900-01		
	.050 Hex Driver Short	HTT		
	.050 Hex Driver Long	HTLT		
	Handpiece Implant Driver Short	DMSIH		
	Handpiece Implant Driver Long	DMLIH		
	Ratchet Implant Driver Short	DTSIH		
	Ratchet Implant Driver Long	DTLIH		
	Thumb Knob for Drivers	ТК		
	Starter Drill Short	DRILLNSS		
	Starter Drill Long	DRILLNSL		
	Drill 2.2mm Short	DRILLN22S		
	Drill 2.2mm Long	DRILLN22L		
	Final Drill for the 3.7mm Implant, short	DRILLN35S		

· · ·	
Surgical	Instruments
Jui yitai	

Size	Part #
Final Drill for the 3.7mm Implant, long	DRILLN35L
Profile Drill for the 3.7mm Implant	DRILLPH37
Final Drill for the 4.1mm Implant, short	DRILL3630S
Final Drill for the 4.1mm Implant, long	DRILL3630L
Profile Drill for the 4.1mm Implant	DRILLPH41
Final Drill for the 4.7mm Implant, short	DRILL4438S
Final Drill for the 4.7mm Implant, long	DRILL4438L
Profile Drill for the 4.7mm Implant	DRILLPH47
Final Drill for the 5.7mm Implant, short	DRILLH53S
Final Drill for the 5.7mm Implant, long	DRILLH53L
Profile Drill for the 5.7mm Implant	DRILLPH57
Final Drill for the 7.0mm Implant, short	DRILL6653S
Profile Drill for the 7.0mm Implant	DRILLPH70
Instrument Box	IHIB
Internal Hex Complete Surgical Kit	IHCSK

### **Internal Hex Drill Taps**

Size	Part #
Tap 3.7mm	TAPH37
Tap 4.1mm	TAPH41
Tap 4.7mm	TAPH47
TAP 5.7mm	TAPH57

	Size	Part #
	Lindemann Bur	LINDR-S21
Ø1.5	Lance Drillt	LD150
	Trephine Bur	TPHB-B6
ĭ	Final Drill for 7.0mm Implant, short	DRILL6653S
	Angled Screw Driver (20Ncm max torque)	SDUGA
	Screw Removal Tool	AFSR

# **Safety Stop Surgical Kit**



	Safety Stop Kit Individual Components And Accessories		
	Size	Part #	
Ø2.0 -	Drill Stop 2.0 x 6mm	DRILNSTP2006	
Ø2.0 -	Drill Stop 2.0 x 8mm	DRILNSTP2008	
Ø2.0 -	Drill Stop 2.0 x 10mm	DRILNSTP2010	
Ø2.0 -	Drill Stop 2.0 x 11.5mm	DRILNSTP2011	
Ø2.0 -	Drill Stop 2.0 x 13mm	DRILNSTP2013	
Ø2.5 -	Drill Stop 2.5 x 6mm	DRILNSTP2506	
<b>⊘</b> 2.5 ∋	Drill Stop 2.5 x 8mm	DRILNSTP2508	
<b>──</b> Ø2.5 ∋	Drill Stop 2.5 x 10mm	DRILNSTP2510	
Ø2.5 =	Drill Stop 2.5 x 11.5mm	DRILNSTP2511	
Ø2.5 -	Drill Stop 2.5 x 13mm	DRILNSTP2513	
Ø3.0 -	Drill Stop 3.0 x 6mm	DRILNSTP3006	
Ø3.0 -	Drill Stop 3.0 x 8mm	DRILNSTP3008	
Ø3.0 -	Drill Stop 3.0 x 10mm	DRILNSTP3010	
Ø3.0 -	Drill Stop 3.0 x 11.5mm	DRILNSTP3011	

	Safety Stop Kit Individual Components And Accessories		
	Size	Part #	
Ø3.0 -	Drill Stop 3.0 x 13mm	DRILNSTP3013	
Ø3.6 -	Drill Stop 3.6 x 6mm	DRILNSTP3606	
Ø3.6 -	Drill Stop 3.6 x 8mm	DRILNSTP3608	
Ø3.6 -	Drill Stop 3.6 x 10mm	DRILNSTP3610	
Ø3.5 -	Drill Stop 3.6 x 11.5mm	DRILNSTP3611	
Ø3.6 -	Drill Stop 3.6 x 13mm	DRILNSTP3613	
Ø4.4 -	Drill Stop 4.4 x 6mm	DRILNSTP4406	
Ø4.4 -	Drill Stop 4.4 x 8mm	DRILNSTP4408	
Ø4.4 -	Drill Stop 4.4 x 10mm	DRILNSTP4410	
Ø4.4 -	Drill Stop 4.4 x 11.5mm	DRILNSTP4411	
Ø4.4 -	Drill Stop 4.4 x 13mm	DRILNSTP4413	
	Drill Stop Drill Extender	TDE	

# **Fully Guided Keyless Kit**

	Part #	
	Fully Guided Keyless Kit	GSKIH
	Master Cylinder for Fully Guided Kit (Cerec M Master Tube Compatible)	TUBE524-504-4L
	Master Cylinder for Fully Guided Kit (Cerec M Master Tube Compatible)	TUBE59-504-45L
	Master Cylinder for Fully Guided Kit (Cerec L Master Tube Compatible)	TUBE524-504-4LNF
	Guided Kit Individual Com And Accessories	ponents
	Size	Part #
Ø2.0 =	Guide Drill 2.0 x 6mm	DRILG2006
Ø2.0 =	Guide Drill 2.0 x 8mm	DRILG2008
Ø2.0 =	Guide Drill 2.0 x 10mm	DRILG2010
Ø2.0	Guide Drill 2.0 x 11.5mm	DRILG2011

	Guided Kit Individual Comp And Accessories	onents
	Size	Part #
Ø2.0 =]	Guide Drill 2.0 x 13mm	DRILG2013
Ø2.5 E	Guide Drill 2.5 x 6mm	DRILG2506
Ø2.5 -	Guide Drill 2.5 x 8mm	DRILG2508
Ø2.5 -	Guide Drill 2.5 x 10mm	DRILG2510
Ø2.5	Guide Drill 2.5 x 11.5mm	DRILG2011
Ø2.5 =	Guide Drill 2.5 x 13mm	DRILG2513
Ø3.0 -	Guide Drill 3.0 x 6mm	DRILG3006
Ø3.0 :	Guide Drill 3.0 x 8mm	DRILG3008
Ø3.0 - ]	Guide Drill 3.0 x 10mm	DRILG3010
Ø3.0 -	Guide Drill 3.0 x 11.5mm	DRILG3011
Ø3.0 -	Guide Drill 3.0 x 13mm	DRILG3013
Ø3.6 -	Guide Drill 3.6 x 6mm	DRILG3606
Ø3.6 -	Guide Drill 3.6 x 8mm	DRILG3608
Ø3.6 =	Guide Drill 3.6 x 10mm	DRILG3610
Ø3.6 =	Guide Drill 3.6 x 11.5mm	DRILG3611
Ø3.6 -	Guide Drill 3.6 x 13mm	DRILG3613
Ø4.4 -	Guide Drill 4.4 x 6mm	DRILG4406
Ø4.4 -	Guide Drill 4.4 x 8mm	DRILG4408
Ø4.4 -	Guide Drill 4.4 x 10mm	DRILG4410
Ø4.4 -	Guide Drill 4.4 x 11.5mm	DRILG4411
Ø4.4 -	Guide Drill 4.4 x 13mm	DRILG4413
	Bone Crest Flattening Drill	DRILGFLAT
	Guided Implant Driver 8.5mm Offset	DTGIH85
	Guided Implant Driver Variable Offset	DTGIHVR
	Guide Fixation Pin (5 Pack)	PIN

### **Hammerless Crestal Sinus Lift Kit**

	Sinus Lift Kit	
	Size	Part #
Vanuo	Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit
	Addition Guide Pin	GPINMAXNP

### **BIO | BlueTouch**

Size	Part #
 Touch-Screen Implant Motor (Adjustable Speed 300 – 40,000 RPM; up to 50 Ncm of torque for placement). Does not include handpiece.	BIO-BT

	Size	Part #		Size
Phan A	Saeshin Strong 20:1 Handpiece External Irrigation (button)	ACL(B)-411		Saeshin Strong 1:1 Straight Handpiece External Irrigation
	Saeshin Traus 20:1 Non-Optic	CRB26XX		
P	Handpiece External Irrigation (button)			Nozzle Clip 1:1 Straight Handpiece
	Anthogyr Impulsion 20:1 Contra Angle (brush type); 50 Ncm tolerance	AHP-85P-I	11	Autoclavable Tubing Set for BIO   BlueTouch
1-	Anthogyr Mont Blanc 20:1 Contra Angle (high efficiency); 80 Ncm	AHP-85MB-X		Autoclavable Pump Section for AE-24
	tolerance			Tubing Set
	Anthogyr Mont Blanc 20:1 Contra Angle with Depth Control (high efficiency); 80 Ncm tolerance	AHP-85MB-CX		Synthetic Premium Handpiece Lubricant
				Lubricant Tool Attachment (necessary for Jubricant)

(necessary for lubricant)

**Sutures** 

All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.

	Size	Part #		Size	Part #
	PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2	version to the state of the st	Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3
	PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2	we 3-0 the 1 second task the	Silk 3-0, FS-2 needle, 12/box	S30FS2
States of the second se	Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3			

# **Master Bone Expander Kit**

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.



ter Bone Expander Kit	
Size	Part #
Master Bone Expander Kit	MBEK



Part #
BEK

#### Lerner Ridge Splitting Chisels

Size	Part #
3mm Chisel with 2mm blade thickness	RSC32
3mm Chisel with .5mm blade thickness	RSC35
6mm Chisel with 2mm blade thickness	RSC62
6mm Chisel with .5mm blade thickness	RSC65

#### **Guided Sinus Lift Kit**

	Guided Sinus Lift Kit	
Littletettet	Size	Part #
	Guided Sinus Lift Kit	Sinus Kit

### **Safety Lateral Sinus Kit**

	Safety Lateral Sinus Kit	
Sading Lations Shour Kill	Size	Part #
	Safety Lateral Sinus Kit	SLSK

#### **Guide Tubes**

Size	Part #
.2 Guide Tube (5 Pack) D=1.2mm DD=2.2mm leight=4mm ip=1mm of total height	TUBE22-12-4L
Guide Tube (5 Pack) .4mm 2.4mm ht=4mm 1mm of total height	TUBE24-14-4L

		Guide Tubes		
Size	Part #		Size	Part #
2.31 Guide Tube (5 Pack) ID=2.31mm OD=3.31mm Height=5mm Lip=1mm of total height	Tube231		2.8 Guide Tube ID=2.8mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-28-5L
2.8 Guide Tube ID=2.8mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-28-5L		3.3 Guide Tube ID=3.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-33-5L
2.9 Guide Tube ID=2.9mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-29-5L		3.5 Guide Tube ID=3.5mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L		3.5 Guide Tube ID=3.5mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L		3.6 Guide Tube ID=3.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-36-5L
2.92 Guide Tube (5 Pack) ID=2.92mm OD=3.92mm Height=5mm Lip=1mm of total height	TUBE292		3.68 Guide Tube (5 Pack) 3.5mm Short/Long, Bio Trilobe ID=3.68mm OD=4.68mm Height=5mm Lip=1mm of total height	TUBE368
3.0 Guide Tube ID=3.0mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-3-5L		3.75 Guide Tube ID=3.75mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-375
3.15 Guide Tube ID=3.15mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-315		3.79 Guide Tube ID=3.79mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-379
3.3 Guide Tube ID=3.3mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-33-5L		3.9 Guide Tube ID=3.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-39-5L
4.1 Guide Tube ID=4.1mm OD=5.0mm Height=4mm Lip=1mm of total height	TUBE5-41-4L		4.45 Guide Tube (5 Pack) 4.3mm Short / Long, Bio Trilobe ID=4.45mm OD=5.45mm	TUBE445

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Zimmer Dental<sup>®</sup>.

		Guide Tubes		
Size	Part #		Size	Part #
4.1 Guide Tube (5 Pack) ID=4.1mm OD=5.1mm Height=3.5mm Lip=1mm of total height	TUBE51-41-35L		4.52 Guide Tube ID=4.52mm OD=5.57mm Height=4mm Lip=1mm of total height	TUBE557-452-4L
4.15 Guide Tube (5 Pack) ID=4.15mm OD=5.15mm Height=5mm Lip=1mm of total height	TUBE415		4.68 Guide Tube ID=4.68mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-468
4.2 Guide Tube ID=4.2mm OD=5.0mm Height=4.5mm Lip=1mm of total height	TUBE5-42-45L		4.7 Guide Tube (5 Pack) ID=4.7mm OD=5.7mm Height=4mm Lip=1mm of total height	TUBE57-47-4L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=4mm Lip=1mm of total height	TUBE525-425-4L		5.0 Guide Tube (5 Pack) ID=5.0mm OD=6.0mm Height=3.5mm Lip=1mm of total height	TUBE6-5-35L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=5mm Lip=1mm of total height	TUBE425		5.0 Guide Tube ID=5.0mm OD=6.2mm Height=4mm Lip=1mm of total height	TUBE62-5-4L
4.25 Guide Tube ID=4.25mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-425		5.02 Guide Tube (5 Pack) ID=5.02mm OD=6.02mm Height=5mm Lip=1mm of total height	TUBE502
4.35 Guide Tube ID=4.35mm OD=5.27mm Height=6mm Lip=1mm of total height	TUBE527-435-6L		5.06 Guide Tube (5 Pack) ID=5.06mm OD=6.02mm Height=3.5mm Lip=1mm of total height	TUBE6-506-35L
5.08 Guide Tube (5 Pack) ID=5.08mm OD=6.08mm Height=5mm Lip=1mm of total height	TUBE508		5.5 Guide Tube ID=5.5mm OD=6.5mm Height=4mm Lip=1mm of total height	TUBE65-55-4L
5.1 Guide Tube ID=5.1mm OD=6.1mm Height=6mm Lip=1mm of total height	TUBE61-51-6L		5.7 Guide Tube (5 Pack) ID=5.7mm OD=6.7mm Height=4mm Lip=1mm of total height	TUBE67-57-4L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L		6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L

		Guide Tubes		
Size	Part #		Size	Part #
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L		6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
5.16 Guide Tube (5 Pack) 5.0mm Short, Bio Trilobe ID=5.16mm OD=6.16mm Height=5mm Lip=1mm of total height	TUBE516		6.2 Guide Tube ID=6.2mm OD=6.9mm Height=4mm Lip=1mm of total height	TUBE69-62-4L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=4mm Lip=1mm of total height	TUBE625-525-4L		6.2 Guide Tube (5 Pack) ID=6.2mm OD=7.2mm Height=3.5mm Lip=1mm of total height	TUBE72-62-35L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=5mm Lip=1mm of total height	TUBE525		6.26 Guide Tube (5 Pack) ID=6.26mm OD=7.21mm Height=3.5mm Lip=1mm of total height	TUBE72-626-35L
5.31 Guide Tube ID=5.31mm OD=6.2mm Height=4.5mm Lip=1mm of total height	TUBE62-531-45L		6.3 Guide Tube ID=6.3mm OD=7.3mm Height=6mm Lip=1mm of total height	TUBE73-63-6L
5.35 Guide Tube (5 Pack) ID=5.35mm OD=6.35mm Height=4mm Lip=1mm of total height	TUBE635-535-4L		6.72 Guide Tube ID=6.72mm OD=7.17mm Height=5mm Lip=1mm of total height	TUBE717-672

#### Terms & Conditions – Important User Information

#### Blue Sky Bio, LLC - TERMS AND CONDITIONS - Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under very limited circumstances.

#### NO RETURNS

Blue Sky Bio does not accept returns. Contact us before attempting to return products. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate. All returned products are subject to inspection prior to any credit/refund authorization. Products returned, that do not meet the exceptions listed above, will incur a restocking fee and may be subject to no credit/refund authorization.

A Return Authorization Number (RA#) must be obtained before returning products.

#### ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

#### SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 11am CST will ship the same day. Orders placed after 11am CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays, as this will affect shipping delivery dates.

#### **INTERNATIONAL BROKER / TAX TERMS**

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

#### CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am - 5pm CST at, 888-446-6724, or email us anytime at, info@blueskybio.com.

#### Important Product Notices For Blue Sky Bio Customers

Blue Sky Bio internal hex dental implants are compatible with Zimmer Dental's Tapered Screw Vent<sup>®</sup> prosthetics, and instrumentation. Blue Sky Bio implants have a patent-pending thread design and should not be inserted with mounts from other companies.

With the exceptions listed above no other companies products should be used with our system.

Federal Law restricts these devices to sale by or on order of a licensed dentist or physician.

Blue Sky Bio is not affiliated with Zimmer Dental.

Tapered Screw Vent<sup>®</sup> is a registered trademark of Zimmer Dental.