


**Your Partner In Digital Dentistry** 

# BIO Conus 12



2020

# Product Catalog BIO Conus 12

Astra OsseoSpeed™ Compatible

January 2020

# **Table of contents**

Implant Packs	2
Practice Building Kits	2,3
Healing Abutments	3
Solid Abutments	4
MPA Custom Abutments	4,5
Straight Abutments	5
KISS Abutments	5
Gold Castable Abutments	5,6
Plastic Castable Abutments	6
Angled Abutments	6,7
Multi Unit Abutments	7
Digital Abutments	8
Intraoral Scanning Parts	8
Overdenture Abutments	8,9
Angled Overdenture Abutments	9
Overdenture Products	10
Analog	10

Implant Level Impression Transfers	10,11
Abutment Screws	11
Surgical Instruments	11-13
Other Surgical Instrumentation	13
Hammerless Crestal Sinus Lift Kit	13
BIO   BlueTouch	14
Suture	14
Master Bone Expander Kit	14
Guided Sinus Lift Kit	15
Safety Lateral Sinus Kit	15
Guide Tubes	15-19
Terms & Conditions	20

### **Implant Packs**

• Packs include Conus 12 implant, 1mm cover screw and a 3mm healing abutment.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.



## **Practice Building Kits**

• Kits include Conus 12 implant, 1mm cover screw, 3mm healing abutment, final abutment/implant level impression coping, and implant level analog.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.





All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Astra®.



### **Healing Abutments**

5.0mm - 5mm height

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.



HDH5005

6.5mm - 5mm height

HDH506505



#### **Solid Abutments**





3.5/4.0 Platform Solid Abutment	
Part #	
ASADHW5052	
NA55WK	
A4WRK	

#### **MPA Custom Abutments**

SMPA Multi Purpose Abutments: Abutments can be used as an implant level impression transfer and the final abutment. Kits include MPA Abutment, Analog and Fixation Screw.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.



3.0 Platform Abutment	
Size	Part #
MPA 3.0 Kit	KMPADH30
MPA 3.0 Abutment	MPADH30
3.0 Implant Level Analog	BDH30



3.5/4.0 Platform Narrow Abutment	
Size	Part #
MPA 3.5/4.0 Regular Kit	KMPADH40
MPA 3.5/4.0 Regular Abutment	MPADH40
3.5/4.0 Implant Level Analog	BDH35






### **Straight Abutments**

Fixation screw included.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra OsseoSpeed™ implant system.





4.5/5.0 Platform Abutment		
111	Size	Part #
	Abutment Conus 12 4.5/5.0mm Straight	ADH5010

#### **KISS Abutments**

KISS Complete Kits include KISS abutment, screw, metal analog. two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

KISS Impression Kits include two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

Fixation screw included.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra OsseoSpeed™ implant system.



3.5/4.0 Regular Platform KISS	
Size	Part #
3.5/4.0 KISS Abutment 1.5mm collar	ASDH35015
Additional Abutment Screw - Conus 12, 3.5mm	ASDH3502

#### **Gold Casatable Abutments**

Gold Casatable Abutments should not be used with high temperture high Palladium content alloys. These are designed for custom abutments using low temperture inlay/onlay gold for cemented crowns. These UCLA abutments are not designed for screw retained crowns. Please see the KISS abutment section for whiteMetal abutment screw retained restorations.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Astra®.



Gold 3.0 Platform Abutment

Size Part #

3.0mm Gold UCLA Abutment
(including screw) AUCLADH30

Additional Abutment Screw, UCLA,
Conus 12, 3.0mm SDH30



3.5/4.0 Platform Abutment

Size Part #


3.5/4.0mm Gold UCLA Abutment AUCLADH35

Additional Abutment Screw, UCLA,
Conus 12, 3.5/4.0mm



4.5/5.0 Platform Abutment	
Size	Part #
4.5/5.0mm Gold UCLA Abutment	AUCLADH50
Additional Abutment Screw, UCLA, Conus 12 4 5/5 0mm	SDH50

#### **Plastic Castable Abutments**





Size	Part #
Plastic 3.0mm UCLA Abutment Non-Engaging (including screw)	PAUCLDHNE30
Additional Abutment Screw, UCLA, Conus 12, 3.0mm	SDH30

Plastic 3.0 Platform Abutment (Non-Engaging)

Size	Part #
Plastic 3.5/4.0mm UCLA Abutment Non-Engaging (including screw)	PAUCLDHNE35
Additional Abutment Screw, UCLA, Conus 12, 3.5/4.0mm	SDH35

Plastic 3.5/4.0 Platform Abutment (Non-Engaging)



Size	Part #
Plastic 3.0mm UCLA Abutment Non-Engaging (including screw)	PAUCLDHNE50
Additional Abutment Screw, UCLA, Conus 12, 3.0mm	SDH50

Plastic 4.5/5.0 Platform Abutment (Non-Engaging)

### **Angled Abutments**

Fixation screw included.

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra Tech Dental® implant system.



3.5/4.0 Platform Angled Abutments	
Size	Part #
Angled Abutment 15 degree	AADH3515



4.5/5.0 Platform Angled Abutments	
Size	Part #
Angled Abutment 15 degree	AADH3515

#### 3.5/4.0 Platform Angled Abutments





4.5/5.0 Platform Angled Abutments	
Size	Part #
Angled Abutment 25 degree	AADH3525

#### **Multi Unit Abutments**

3.5/4.0 Platform	
Size	Part #
Abutment Multi Unit Straight, 3.5/4.0 platform, Conus 12, 1.5mm cuff - Requires .080 Hex Driver	AMUSDH3515
Abutment Multi Unit Straight, 3.5/4.0 platform, Conus 12, 3.0mm cuff - Requires .080 Hex Driver	AMUSDH3530
Abutment Multi Unit Straight, 3.5/4.0 platform, Conus 12, 5.0mm cuff - Requires .080 Hex Driver	AMUSDH3550
Angled Abutment Multi Unit, 3.5/4.0 platform, 17 degree, Conus 12, includes abutment screw	AAMUDH3517
Angled Abutment Multi Unit, 3.5/4.0 platform, 30 degree, Conus 12, includes abutment screw	AAMUDH3530



	Com	ponents		
Size	Part #		Size	Part #
Comfort Cap for all Multi Unit Abutments	СМИ		MMUTC Casting Sleeve	PWAXMUTC
Combination Impression Transfer, Temporary Abutment and Castable Sleeve, Plastic - Requires screw SMUS	MMUP		Analog, Multi Unit	BMU
or SMUL Combination Impression Transfer/ Temporary Abutment - Titanium Alloy - Requires screw SMUS or SMUL	MMUT		Angled Multi Unit Abutment Delivery Aid  Screw Multi Unit Impression Transfer, long	DSMU SMUL
Multiunit Coping, Temporary Abutment Titanium Alloy - Requires screw SMUS or SMUL	MMUTC	<u>_</u>	Screw, Multi Unit, Castable Abutments, short	SMUS
Multi Unit Titanium Coping Thimble	MUT		.080 Hex Driver for straight Multi Unit Abutments	HFD
MMUTC Spacer Sleeve	PSPACEMUTC			

### **Digital Abutment**



Narrow Platform	
Size	Part #
Digital Abutment for abutment and crowns for zirconia and emax (includes fixation screw)	ABDH30
Additional Abutment Screw	SDH30



Regular Platform	
Size	Part #
Digital Abutment for abutment and crowns for zirconia and emax (includes fixation screw)	ABDH35
Additional Abutment Screw	SDH35



Regular Platform		
Size	Part #	
Digital Abutment for abutment and crowns for zirconia and emax (includes fixation screw)	ABDH50	
Additional Abutment Screw	SDH50	

## **Intraoral Scanning Parts**

#### Master Overdenture Restorative Kits Regular - Chairside

Size	Part #
Rosetta Scanpost Extender for Regular Digital Abutments	OPTIBASER
Regular Ball ScanSnap (required for Exocad or 3Shape)	BALL6
Rosetta Scanpost Extender for Mini Digital Abutments	OPTIBASEN

#### Master Overdenture Restorative Kits Regular - Laboratory

Size	Part #
Mini Ball ScanSnap (required for Exocad / 3Shape)	BALL6MI
ScanFlag Scanpost for Regular Digital Abutments	OPTICSTB
ScanFlag Scanpost for Mini Digital Abutments	OPTICSTBMI

#### **Overdenture Abutments**

Super Snap Overdenture Restorative Kit – Chairside includes: overdenture abutment, retention inserts and block out ring Super Snap Overdenture Restorative Kit – Laboratory includes: overdenture abutment, retention inserts, block out ring, analog and impression transfer

Super Snap Overdenture Restorative Kits 3.0 platform - ChairsideRequires No Metal Housing

Size	Part #
Super Insert Restorative Kits - Chairside, regular x 3.0mm height	SIRKDH3003

Super Insert Restorative Kit
- Chairside, regular x 5.0mm height

Super Snap Overdenture Restorative Kits 3.0 platform - LaboratoryRequires No Metal Housing



Size	Part #
Super Insert Restorative Kits - Laboratory, regular x 3.0mm height	SIRKDH3003-L



Super Insert Restorative Kits
- Laboratory, regular x 5.0mm height

SIRKDH3005-L

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Astra<sup>®</sup>.

#### **Super Snap Overdenture Restorative Kits 3.5/4.0 platform** - ChairsideRequires No Metal Housing

Size	Part #
Super Insert Restorative Kits - Chairside, regular x 3.0mm height	SIRKDH3503



Super Insert Restorative Kit

- Chairside, regular x 5.0mm height

SIRKDH3505

#### **Super Snap Overdenture Restorative Kits 3.5/4.0 platform** LaboratoryRequires No Metal Housing





Super Insert Restorative Kits

- Laboratory, regular x 5.0mm height

SIRKDH3505-L

#### **Super Snap Overdenture Restorative Kits 4.5/5.0 platform** - ChairsideRequires No Metal Housing

Size	Part #
Super Insert Restorative Kits - Chairside, regular x 3.0mm height	SIRKDH5003
Constitution of Destauration 1/4	



Super Insert Restorative Kit

- Chairside, regular x 5.0mm height

SIRKDH5005

#### **Super Snap Overdenture Restorative Kits 4.5/5.0 platform** - LaboratoryRequires No Metal Housing

Size	Part #
Super Insert Restorative Kits - Laboratory, regular x 3.0mm height	SIRKDH5003-L



Super Insert Restorative Kits

- Laboratory, regular x 5.0mm height

SIRKDH5005-L

### **Angled Overdentures**

#### Angled Super Overdenture Restorative Kit 3.5/4.0 Platform - Chairside Requires No Metal Housing

	•	
	Size	Part #
	Angled Super Overdenture Restorative Kit Narrow Platform, 17 degrees - Chairside	ASRKC3517
	Angled Super Overdenture Restorative Kit Narrow Platform, 30 degrees	ASDV02520



- Chairside

- Chairside

ASRKC3530

#### Angled Super Overdenture Restorative Kit 3.5/4.0 Platform - Laboratory Requires No Metal Housing

		9
	Size	Part #
	Angled Super Overdenture Restorative Kit Narrow Platform, 17 degrees - Laboratory	ASRKC3517-L
4	Angled Super Overdenture Restorative Kit Narrow Platform, 30 degrees - Laboratory	ASRKC3530-L




	I latioiiii	- Chairside Requires No Metai Ti	odsing
		Size	Part #
0		Angled Super Overdenture Restorative Kit Regular Platform, 17 degrees - Chairside	ASRKSQR17
		Angled Super Overdenture Restorative Kit Regular Platform, 30 degrees	ASRKSQR30

ASRKSQR30

#### **Angled Super Overdenture Restorative Kit 4.5/5.0 Platform - Laboratory Requires No Metal Housing**

Size	Part #
Angled Super Overdenture Restorative Kit Regular Platform, 17 degrees - Laboratory	ASRKSQR17-L



Angled Super Overdenture Restorative Kit RegularPlatform, 30 degrees

ASRKSQR30-L

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Astra®.

### **Overdenture Products**





### **Analogs**

Size	Part	#	Size		Part #
3.0 Implant Level Analog	BDH3	0	3.5/4.0	) Implant Level Analog	BDH35
	Size		Part #	I	
	4.5/5.0 Implant Le	evel Analog	BDH50		

# **Implant Level Impression Transfers**








#### **Abutment Screws**



## **Surgical Instruments**

BIO | Conus 12 implants, prosthetics and instrumentation are fully compatible with Astra OsseoSpeed™ implant system.






All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Astra®.

	Surgical Instruments	
	Size	Part #
	Final Drill for the 4.5mm Implant, short	DRILLD42S
	Final Drill for the 4.5mm Implant, long	DRILLD42L
	Profile Drill for the 4.5mm Implant	DRILLPDH45
	Final Drill for the 5.0mm Implant, short	DRILLD46S
	Final Drill for the 5.0mm Implant, long	DRILLD46L
	Profile Drill for the 5.0mm Implant	DRILLPDH50
	Instrument Box	C12IB
A COMMANDA OF THE COMMAND OF THE COM	Conus 12 Complete Surgical Kit	C12CSK

# **Other Surgical Instrumentation**



### **Hammerless Crestal Sinus Lift Kit**



### **BIO | BlueTouch**



Size

BIO-BT

Part #



Size	Part #
Saeshin Strong 20:1 Handpiece External Irrigation (button)	ACL(B)-411
Saeshin Traus 20:1 Non-Optic	CRB26XX
Handpiece External Irrigation (button)	CNB20XX
Anthogyr Impulsion 20:1 Contra Angle (brush type); 50 Ncm tolerance	AHP-85P-I
Anthogyr Mont Blanc 20:1 Contra Angle (high efficiency); 80 Ncm tolerance	AHP-85MB-X
Anthogyr Mont Blanc 20:1 Contra Angle with Depth Control (high efficiency); 80 Ncm tolerance	AHP-85MB-C



### **Sutures**

All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.





## **Master Bone Expander Kit**

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.




Size	Part #
Master Bone Expander Kit	MBEK

r Kit
Part #
BEK

Lerner Ridge Splitting Chisels					
1111	Size	Part #			
	3mm Chisel with 2mm blade thickness	RSC32			
	3mm Chisel with .5mm blade thickness	RSC35			
1111	6mm Chisel with 2mm blade thickness	RSC62			
1111	6mm Chisel with .5mm blade thickness	RSC65			

All prices are available in the product section of our website. Volume discounts are available for some items. Blue Sky Bio is not affiliated with Astra®.

## **Guided Sinus Lift Kit**



# **Safety Lateral Sinus Kit**

F	Safety Lateral Sinus Kit	
Select Lateral History	Size	Part #
	Safety Lateral Sinus Kit	SLSK

Guide Tubes				
Size	Part #	Size	Part #	
1.2 Guide Tube (5 Pack) ID=1.2mm OD=2.2mm Height=4mm Lip=1mm of total height	TUBE22-12-4L	1.9 Guide Tube (5 Pack) ID=1.9mm OD=2.9mm Height=4mm Lip=1mm of total height	TUBE29-19-4	
1.4 Guide Tube (5 Pack) ID=1.4mm OD=2.4mm Height=4mm Lip=1mm of total height	TUBE24-14-4L	2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=4mm Lip=1mm of total height	TUBE3-2-4l	
1.55 Guide Tube (5 Pack) ID=1.55mm OD=2.55mm Height=7mm		2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=5mm Lip=1mm of total height	Tube20	
1.6 Guide Tube (5 Pack) ID=1.6mm OD=2.6mm Height=4mm Lip=1mm of total height	TUBE26-16-4L	2.1 Guide Tube (5 Pack) ID=2.1mm OD=3.1mm Height=4mm Lip=1mm of total height	TUBE31-21-4	
1.7 Guide Tube (5 Pack) ID=1.7mm OD=2.7mm Height=4mm Lip=1mm of total height	TUBE27-17-4L	2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-21-5L	

Size	Part #	Size	Part #
2.1 Guide Tube ID=2.1mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-21-5L	2.4 Guide Tube (5 Pack) ID=2.4mm OD=3.4mm Height=4mm Lip=1mm of total height	TUBE34-24-4L
2.1 Guide Tube ID=2.1mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-21-5L	2.4 Guide Tube ID=2.4mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-24-5L
2.11 Guide Tube (5 Pack) ID=2.11mm OD=3.11mm Height=5mm Lip=1mm of total height	Tube211	2.4 Guide Tube ID=2.4mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-24-5L
2.2 Guide Tube (5 Pack) ID=2.2mm OD=3.2mm Height=5mm Lip=1mm of total height	Tube22	2.4 Guide Tube ID=2.4mm OD=5.0mm Height=5mm Lip=1mm of total heigh	TUBE5-24-5L
2.3 Guide Tube ID=2.3mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-23-5L	2.6 Guide Tube ID=2.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-26-5L
2.3 Guide Tube ID=2.3mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-23-5L	2.72 Guide Tube ID=2.72mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-272
2.3 Guide Tube ID=2.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-23-5L	2.8 Guide Tube (5 Pack) ID=2.8mm OD=3.8mm Height=5mm Lip=1mm of total height	TUBE280
2.31 Guide Tube (5 Pack) ID=2.31mm OD=3.31mm Height=5mm Lip=1mm of total height	Tube231	2.8 Guide Tube ID=2.8mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-28-5L
2.8 Guide Tube ID=2.8mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-28-5L	3.3 Guide Tube ID=3.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-33-5L
2.9 Guide Tube ID=2.9mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-29-5L	3.5 Guide Tube ID=3.5mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-35-5L

Size	Part #	Size	Part #
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.5 Guide Tube ID=3.5mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-35-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L	3.6 Guide Tube ID=3.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-36-5L
2.92 Guide Tube (5 Pack) ID=2.92mm OD=3.92mm Height=5mm Lip=1mm of total height	TUBE292	3.68 Guide Tube (5 Pack) 3.5mm Short/Long, Bio Trilobe ID=3.68mm OD=4.68mm Height=5mm Lip=1mm of total height	TUBE368
3.0 Guide Tube ID=3.0mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-3-5L	3.75 Guide Tube ID=3.75mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-375
3.15 Guide Tube ID=3.15mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-315	3.79 Guide Tube ID=3.79mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-379
3.3 Guide Tube ID=3.3mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-33-5L	3.9 Guide Tube ID=3.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-39-5L
4.1 Guide Tube ID=4.1mm OD=5.0mm Height=4mm Lip=1mm of total height	TUBE5-41-4L	4.45 Guide Tube (5 Pack) 4.3mm Short / Long, Bio Trilobe ID=4.45mm OD=5.45mm	TUBE445
4.1 Guide Tube (5 Pack) ID=4.1mm OD=5.1mm Height=3.5mm Lip=1mm of total height	TUBE51-41-35L	4.52 Guide Tube ID=4.52mm OD=5.57mm Height=4mm Lip=1mm of total height	TUBE557-452-4L
4.15 Guide Tube (5 Pack) ID=4.15mm OD=5.15mm Height=5mm Lip=1mm of total height	TUBE415	4.68 Guide Tube ID=4.68mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-468
4.2 Guide Tube ID=4.2mm OD=5.0mm Height=4.5mm Lip=1mm of total height	TUBE5-42-45L	4.7 Guide Tube (5 Pack) ID=4.7mm OD=5.7mm Height=4mm Lip=1mm of total height	TUBE57-47-4L

Size	Part #	Size	Part #
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=4mm Lip=1mm of total height	TUBE525-425-4L	5.0 Guide Tube (5 Pack) ID=5.0mm OD=6.0mm Height=3.5mm Lip=1mm of total height	TUBE6-5-35L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=5mm Lip=1mm of total height	TUBE425	5.0 Guide Tube ID=5.0mm OD=6.2mm Height=4mm Lip=1mm of total height	TUBE62-5-4L
4.25 Guide Tube ID=4.25mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-425	5.02 Guide Tube (5 Pack) ID=5.02mm OD=6.02mm Height=5mm Lip=1mm of total height	TUBE502
4.35 Guide Tube ID=4.35mm OD=5.27mm Height=6mm Lip=1mm of total height	TUBE527-435-6L	5.06 Guide Tube (5 Pack) ID=5.06mm OD=6.02mm Height=3.5mm Lip=1mm of total height	TUBE6-506-35L
5.08 Guide Tube (5 Pack) ID=5.08mm OD=6.08mm Height=5mm Lip=1mm of total height	TUBE508	5.5 Guide Tube ID=5.5mm OD=6.5mm Height=4mm Lip=1mm of total height	TUBE65-55-4L
5.1 Guide Tube ID=5.1mm OD=6.1mm Height=6mm Lip=1mm of total height	TUBE61-51-6L	5.7 Guide Tube (5 Pack) ID=5.7mm OD=6.7mm Height=4mm Lip=1mm of total height	TUBE67-57-4L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L	6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
5.16 Guide Tube (5 Pack) 5.0mm Short, Bio Trilobe ID=5.16mm OD=6.16mm Height=5mm Lip=1mm of total height	TUBE516	6.2 Guide Tube ID=6.2mm OD=6.9mm Height=4mm Lip=1mm of total height	TUBE69-62-4L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=4mm Lip=1mm of total height	TUBE625-525-4L	6.2 Guide Tube (5 Pack) ID=6.2mm OD=7.2mm Height=3.5mm Lip=1mm of total height	TUBE72-62-35L

Size	Part #	Size	Part #
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=5mm Lip=1mm of total height	TUBE525	6.26 Guide Tube (5 Pack) ID=6.26mm OD=7.21mm Height=3.5mm Lip=1mm of total height	TUBE72-626-35L
5.31 Guide Tube ID=5.31mm OD=6.2mm Height=4.5mm Lip=1mm of total height	TUBE62-531-45L	6.3 Guide Tube ID=6.3mm OD=7.3mm Height=6mm Lip=1mm of total height	TUBE73-63-6L
5.35 Guide Tube (5 Pack) ID=5.35mm OD=6.35mm Height=4mm Lip=1mm of total height	TUBE635-535-4L	6.72 Guide Tube ID=6.72mm OD=7.17mm Height=5mm Lip=1mm of total height	TUBE717-672

### **Terms & Conditions – Important User Information**

#### Blue Sky Bio, LLC - TERMS AND CONDITIONS - Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

#### **NO RETURNS**

Blue Sky Bio does not accept returns. Contact us before attempting to return products. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate. All returned products are subject to inspection prior to any credit/refund authorization. Products returned, that do not meet the exceptions listed above, will incur a restocking fee and may be subject to no credit/refund authorization.

A Return Authorization Number (RA#) must be obtained before returning products.

#### ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

#### SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 11am CST will ship the same day. Orders placed after 11am CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays, as this will affect shipping delivery dates.

#### **INTERNATIONAL BROKER / TAX TERMS**

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

#### CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am - 5pm CST at, 888-446-6724, or email us anytime at, info@blueskybio.com.

#### Important Product Notices For Blue Sky Bio Customers

Blue Sky Bio dental implants are compatible with prosthetics, and instrumentation.

Blue Sky Bio implants have a patent-pending thread design and should not be inserted with mounts from other companies.

With the exceptions listed above no other companies products should be used with our system.

Federal Law restricts these devices to sale by or on order of a licensed dentist or physician.