


# Compatibility Innovation Value

BIO | Max



# 2017

## Product Catalog

BIO | Max

Nobel Active® Compatible

January 2017



Blue Sky Bio



## Table of contents

**Implants**

**2**

**Short Implants**

**2,3**

**3.0mm Implants**

**3**

**Healing Abutments**

**3,4**

**MPA Custom Abutments**

**4**

**KISS Abutments**

**4,5**

**UCLA Gold and Plastic Castable Abutments**

**5**

**Angled Abutments**

**6**

**Multi Unit Abutments**

**6**

**Temporary Abutment**

**7**

**CAD/CAM Tibase**

**7**

**Overdentures - Straight and Angled**

**7,8,9**

**Analogs**

**9**

**Impression Transfers**

**9**



<b>Surgical Instruments</b>	<b>10,11</b>
<b>BIO   BlueTouch</b>	<b>11</b>
<b>Suture</b>	<b>12</b>
<b>Master Bone Expander Kit</b>	<b>12</b>
<b>Sinus Lift Kit</b>	<b>12</b>
<b>Safety Lateral Sinus Kit</b>	<b>13</b>
<b>BIO   Max Safety Stop Surgical Kit</b>	<b>13</b>
<b>Fully Guided Keyless Kit</b>	<b>13</b>
<b>Guided Tubes</b>	<b>14-18</b>
<b>Terms and Conditions</b>	<b>19</b>


## Implants

Packs include BIO | Max implant and 1mm cover screw.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.


BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

**3.5mm Ø NP platform switched**




Size	Part #
8mm	IJH3508
10mm	IJH3510
11.5mm	IJH3511
13mm	IJH3513
16mm	IJH3516

**4.3mm Ø NP platform switched**




Size	Part #
8mm	IJH4308
10mm	IJH4310
11.5mm	IJH4311
13mm	IJH4313
16mm	IJH4316

**5.0mm Ø NP platform switched**



Size	Part #
8mm	IJH5008
10mm	IJH5010
11.5mm	IJH5011
13mm	IJH5013
16mm	IJH5016

**7.0mm Ø NP platform switched**



Size	Part #
8mm	IJH7008
10mm	IJH7010
11.5mm	IJH7011
13mm	IJH7013
16mm	IJH7016


## Short Implants

Kits include BIO | Max implant and 1mm cover screw.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.


BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

**4.3mm Ø NP platform switched**



Size	Part #
6mm	IJH4306

**5.0mm Ø NP platform switched**



Size	Part #
6mm	IJH5006

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

### 7.0mm Ø NP platform switched



Size	Part #
6mm	IJH7006

## 3.0mm Implants

Kits include BIO | Max implant and 1mm cover screw.

Requires DRILLN26S or DRILLN26L.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

### 3.5mm Ø NP platform switched



Size	Part #
10mm	IJH3510
11.5mm	IJH3511
13mm	IJH3513

## Healing Abutments

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

### 3.5 Healing Abutments

	Size	Part #
	3mm height	HJH3503
	5mm height	HJH3505

### 4.3 Healing Abutments

	Size	Part #
	3mm height	HJH4303
	5mm height	HJH4305

### 4.7 Healing Abutments

	Size	Part #
	3mm height	HJH4703
	5mm height	HJH4705

All prices are available in the product section of our website. Volume discounts are available for some items.


Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.



6.0 Healing Abutments		
	Size	Part #
	3mm height	HJH6003
	5mm height	HJH6005




Cover Screw		
	Size	Part #
	Cover Screw	CJH




## MPA Custom Abutments




MPA Multi Purpose Abutments: Abutments can be used as an implant level impression transfer and the final abutment. Kits include MPA Abutment, Analog and Fixation Screw.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

4.1 Abutment NP platform		
	Size	Part #
	MPA Kit	KMPAJH41
	MPA Abutment	MPAJH41
	Implant level Analog	BIJH

4.7 Abutment NP platform		
	Size	Part #
	MPA Kit	KMPAJH47
	MPA Abutment	MPAJH47
	Implant level Analog	BIJH

6.0 Abutment NP platform		
	Size	Part #
	MPA Kit	KMPAJH60
	MPA Abutment	MPAJH60
	Implant level Analog	BIJH

## KISS Abutments

Fixation screw included

KISS Complete Kits include KISS abutment, screw, metal analog. two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

KISS Impression Kits include two resin analogs, impression transfer, engaging and non-engaging waxing sleeves and temporary sleeve.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.



## NP Platform KISS



Size	Part #
KISS Complete Kit 1.5mm collar	KASJHR
KISS Abutment 1.5mm collar	ASJHR015
KISS Abutment Analog	BMKISR
KISS Impression Kit	KIKR

## NP Platform Wide KISS



Size	Part #
KISS Wide Complete Kit 1.5mm collar	KASJHW
KISS Abutment 1.5mm collar, Wide	ASJHW015
KISS Abutment Analog, Wide	BMKISW
KISS Impression Kit Wide	KIKW

## UCLA Gold Abutments

Fixation screw included

Please see the KISS abutment section for whiteMetal abutment screw retained restorations.

BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

### Abutment NP platform



Size	Part #
Gold UCLA Abutment	AUCLAJHG
Additional Abutment Screw, UCLA, BIO   Max	SRJH

### Abutment NP platform



Size	Part #
Gold Non-Engaging UCLA Abutment	AUCLAJHNG
Additional Abutment Screw, UCLA, BIO   Max	SRJH

## Plastic Castable Abutments

### Abutment NP platform



Size	Part #
Plastic UCLA Abutment (including screw)	PAUCLJHNP
Additional Abutment Screw, UCLA, BIO   Max	SRJH

### Abutment NP platform



Size	Part #
Plastic Non-Engaging UCLA Abutment (including screw)	PAUCLJHNP
Additional Abutment Screw, UCLA, BIO   Max	SRJH

### Abutment NP platform



Size	Part #
MPA Plastic Kit Castable Non - Engaging (including abutment, analog and screw)	KPMPAJHNE41
MPA Plastic Abutment Castable Non - Engaging 4.1mm	PMPAJHNE41
Implant level Analog	BIJH

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Angled Abutments





Fixation screw included.









BIO | Max implants are prosthetically compatible with NobelActive® NP platform prosthetic abutments\*\*.

BIO | Max Prosthetics are compatible with NobelActive® NP implant platform.

Angled Abutments		
	Size	Part #
	Angled Abutment 15 degree	AAJH15
	Angled Abutment 30 degree	AAJH30

## Multi Unit Abutments

Multi Unit Abutments					
	Size	Part #		Size	Part #
	Abutment Multi Unit Straight, Max, 3.5mm x 1.5mm cuff - <b>Requires .080 Hex Driver</b>	AMUSJH15		Angled Abutment Multi Unit, Max, 3.5mm 17 degree	AAMUJH17
	Abutment Multi Unit Straight, Max, 3.5mm x 3.0mm cuff - <b>Requires .080 Hex Driver</b>	AMUSJH30		Angled Abutment Multi Unit, Max, 3.5mm 30 degree	AAMUJH30

Components					
	Size	Part #		Size	Part #
	Comfort Cap for all Multi Unit Abutments	CMU		Angled Multi Unit Abutment Delivery Aid	DSMU
	Combination Impression Transfer, Temporary Abutment and Castable Sleeve, Plastic - <b>Requires screw SMUS or SMUL</b>	MMUP		Screw Multi Unit Impression Transfer, long	SMUL
	Combination Impression Transfer/Temporary Abutment - Titanium Alloy - <b>Requires screw SMHS or SMUL</b>	MMUT		Screw, Multi Unit, Castable Abutments, short	SMUS
	Analog, Multi Unit	BMU		.080 Hex Driver for straight Multi Unit Abutments	HFD

All prices are available in the product section of our website. Volume discounts are available for some items.



Blue Sky Bio is not affiliated with Nobel Biocare.



Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.









## Temporary Abutment

	Size	Part #
	Titanium Base Abutment for temporary abutment, Narrow Platform (includes abutment screw)	ABJHNP
	Titanium Base Abutment for temporary abutment 1.8mm collar, Narrow Platform (includes abutment screw)	ABJHNP18

	Size	Part #
	Titanium Base Abutment for temporary abutment, Mini, Narrow Platform (includes abutment screw)	ABJHNPMI
	Titanium Base Abutment for temporary abutment - non engaging, Narrow Platform (includes abutment screw)	ABJHNPN

## CAD/CAM Tibase

CAD/CAM Titanium Base		
	Size	Part #
	Titanium Base Abutment for CAD/CAM abutment and crowns for zirconia and emax (includes fixation screw)	ABJHNP
	Titanium Base Abutment Mini for CAD/CAM abutment and crowns for zirconia and emax (includes fixation screw)	ABJHNPMI
	Titanium Base Abutment for temporary abutment - non engaging, Narrow Platform (includes abutment screw)	ABJHNP

Angled CAD/CAM Titanium Base		
	Size	Part #
	Titanium Angled Base Abutment for temporary abutment, Narrow Platform (includes abutment screw)	AABJHNP
	Titanium Angled Base Abutment for temporary abutment - non engaging, Narrow Platform (includes abutment screw)	AABJHNPN
	Titanium Angled Base Abutment for temporary abutment - without hex, Narrow Platform (includes abutment screw)	AABJHNPNC

## Overdentures

### Master Overdenture Restorative Kits Regular - Chairside

	Size	Part #
	Master Overdenture Restorative Kit - Chairside, regular x 3.0mm height	MORKJH30
	Master Overdenture Restorative Kit - Chairside, regular x 5.0mm height	MORKJH50

### Master Overdenture Restorative Kits Regular - Laboratory

	Size	Part #
	Master Overdenture Restorative Kit - Laboratory, regular x 3.0mm height	MORKJH30-L
	Master Overdenture Restorative Kit - Laboratory, regular x 5.0mm height	MORKJH30-L


All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Super Snap Overdenture Restorative Kits - Chairside Requires No Metal Housing

	Size	Part #
	Super Insert Restorative Kits - Chairside, 3mm	SIRKJH30
	Super Insert Restorative Kits - Chairside, 5mm	SIRKJH50

## Super Snap Overdenture Restorative Kits - Laboratory Requires No Metal Housing

	Size	Part #
	Super Insert Restorative Kits - Laboratory, 3mm	SIRKJH30-L
	Super Insert Restorative Kits - Laboratory, 5mm	SIRKJH50-L

## Angled Overdentures


### Angled Overdenture Restorative Kits Regular - Chairside

	Size	Part #
	Angled Overdenture Restorative Kit , 17 degrees - Chairside	ARKJH17
	Angled Overdenture Restorative Kit, 30 degrees - Chairside	ARKJH30



### Angled Overdenture Restorative Kits Regular - Laboratory

	Size	Part #
	Angled Overdenture Restorative Kit, 17 degrees - Laboratory	ARKJH17-L
	Angled Overdenture Restorative Kit , 30 degrees - Laboratory	ARKJH30-L

### Angled Super Snap Overdenture Restorative Kits - Chairside Requires No Metal Housing

	Size	Part #
	Angled Super Overdenture Restorative Kit , 17 degrees - Chairside	ASRKJH17
	Angled Super Overdenture Restorative Kit , 17 degrees - Chairside	ASRKJH30

### Angled Super Snap Overdenture Restorative Kits - Laboratory Requires No Metal Housing

	Size	Part #
	Angled Super Overdenture Restorative Kit , 17 degrees - Laboratory	ASRKJH17-L
	Angled Super Overdenture Restorative Kit , 17 degrees - Laboratory	ASRKJH30-L








All prices are available in the product section of our website. Volume discounts are available for some items.









Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Overdenture Products


	Size	Part #
	Complete Flat Snap Restorative Kit	NOPK
	Analog	FSA
	Replacement Inserts Strong Retention (4 Pack)	FSCS
	Replacement Inserts Medium / Light Retention (4 Pack)	FSCR
	Resin Impression transfer (2 Pack) and FSCB (4 Pack)	FSI Pack
	Super Insert – Strong retention, requires no housing (includes 2 inserts and 2 block out rings) <b>Extreme long lasting retention on parallel implants.</b>	
	Super Insert – Medium retention, requires no housing (includes 2 inserts and 2 block out rings) <b>Extreme long lasting retention on parallel implants.</b>	


	Size	Part #
	Processing Inserts No Retention (4 Pack)	FSCB
	Resin Impression transfer (2 Pack)	FSI
	Titanium Housing	HGFS
	.050 Hex Driver Short	HTT
	.048 Hex Driver Short	SDS048
	Insert Placement/Removal Instrument	TFSU
	Super Insert – Light retention, requires no housing (includes 2 inserts and 2 block out rings) <b>Extreme long lasting retention on parallel implants.</b>	
	Insert Placement/Removal Instrument	

## Analogs

	Size	Part #
	Implant Level Analog	BIJH

## Impression Transfers

	Size	Part #
	Closed Tray Impression Transfer	MIJH

	Size	Part #
	Open Tray Impression Transfer	MIJH - OT




















All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Surgical Instruments

Surgical Instruments	
Size	Part #
	Adapter for NB Ratchet NA
	Torque/Ratchet TR-900-01
	Screwdriver short SDS048
	Screwdriver long HNLT
	Handpiece Implant Driver Short DMJHS
	Handpiece Implant Driver Long DMJFL
	Ratchet Implant Driver Short DTJHS
	Ratchet Implant Driver Long DTJHL
	Starter Drill Short DRILLNSS
	Starter Drill Long DRILLNSL
	Drill 2.2mm Short DRILLN22S
	Drill 2.2mm Long DRILLN22L
	Implant Drill Short 3.5mm DRILLN35S
	Implant Drill Long 3.5mm DRILLN35L
	Implant Drill Short 4.3mm DRILLN43S
	Implant Drill Long 4.3mm DRILLN43L
	Implant Drill Short 5.0mm DRILLN50S
	Implant Drill Long 5.0mm DRILLN50L
	Implant Drill Short 7.0mm DRILL6653S

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Surgical Instruments



Size	Part #
Instrument Box	MIB

BIO | Max Complete Surgical Kit      MCSK

## Other Surgical Instrumentation







Size	Part #
Lindemann Bur	LINDR-S21

Trepine Bur	TPHB-B6
-------------	---------

## BIO | BlueTouch



Size	Part #
Touch-Screen Implant Motor (Adjustable Speed 300 – 40,000 RPM; up to 50 Ncm of torque for placement). Does not include handpiece.	BIO-BT

Size	Part #
 Saeshin Strong 20:1 Handpiece External Irrigation (button)	ACL(B)-411
 Saeshin Traus 20:1 Non-Optic Handpiece External Irrigation (button)	CRB26XX
 Anthogyr Impulsion 20:1 Contra Angle (brush type); 50 Ncm tolerance	AHP-85P-I
 Anthogyr Mont Blanc 20:1 Contra Angle (high efficiency); 80 Ncm tolerance	AHP-85MB-X
 Anthogyr Mont Blanc 20:1 Contra Angle with Depth Control (high efficiency); 80 Ncm tolerance	AHP-85MB-CX

Size	Part #
 Saeshin Strong 1:1 Straight Handpiece External Irrigation	AT-11
Nozzle Clip 1:1 Straight Handpiece	ST-Nozzle
 Autoclavable Tubing Set for BIO   BlueTouch	AE-24
 Autoclavable Pump Section for AE-24 Tubing Set	AE-24-PUMP
 Synthetic Premium Handpiece Lubricant	LU1011HS
 Lubricant Tool Attachment (necessary for lubricant)	LT1013HS

All prices are available in the product section of our website. Volume discounts are available for some items.






Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.



## Sutures

All sutures are 60cm in length with reverse cutting needles and packaged 12 per box. Sutures sales are currently limited to U.S. customers.


Size	Part #	Size	Part #
 PGA Fast Resorb 4-0, FS-2 needle, 12/box	PGF40FS2	 Nylon (unbraided black) 5-0, P3 needle, 12/box	N50P3
 PGA 4-0 violet, FS-2 needle, 12/box	PG40FS2	 Silk 3-0, FS-2 needle, 12/box	S30FS2
 Polypropylene (unbraided black) 5-0, P3 needle, 12/box	P50P3		

## Master Bone Expander Kit

Master Bone Expander Kit includes the Bone Expander Kit and four Lerner Ridge Splitting Chisels.


Master Bone Expander Kit		Bone Expander Kit	
Size	Part #	Size	Part #
 Master Bone Expander Kit	MBEK	 Bone Expander Kit	BEK

Lerner Ridge Splitting Chisels	
Size	Part #
 3mm Chisel with 2mm blade thickness	RSC32
3mm Chisel with .5mm blade thickness	RSC35
6mm Chisel with 2mm blade thickness	RSC62
6mm Chisel with .5mm blade thickness	RSC65

## Sinus Lift Kit

Complete Sinus Lift Kit includes: Sinus Probe, Bone Carrier, Bone Condenser, 2mm Sinus Drill, Sinus Reamers (2.4mm, 2.8mm, 3.2mm, 3.6mm), Stoppers (2mm – 11mm), Stopper Stands, and instrument box.

Sinus Lift Kit	
Size	Part #
 Hammerless Atraumatic Crestal Sinus Lift Kit	Sinus Kit

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Safety Lateral Sinus Kit

Safety Lateral Sinus Kit	
Size	Part #
Safety Lateral Sinus Kit	SLSK



## BIO | Max Safety Stop Surgical Kit

The Safety Stop Surgical Kit includes the following: 25 drills with integrated stops, 4 implant drivers, 2 .048 screw drivers, 3 paralleling pins, drill extender and ratchet adapter.

Sinus Lift Kit	
Size	Part #
Safety Stop Surgical Kit	SSSMNP



## Fully Guided Keyless Kit

The Guided Surgical Kit includes the following: 25 drills, 2 implant drivers, 2 .048 screw drivers, 3 paralleling pins and Thumb Knob. Handpiece implant drivers are not included.

Surgical Guides that are manufactured by 3rd party vendors with a Nobel Master Tubes can be used with the Fully Guided Keyless Kit. It will allow for implant placement via the guide up to and including a 5.0mm implant.

Sinus Lift Kit	
Size	Part #
Fully Guided Keyless Kit	SSSMNP



All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.

\*\* BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

## Guide Tubes

### Guide Tubes

Size	Part #	Size	Part #
1.2 Guide Tube (5 Pack) ID=1.2mm OD=2.2mm Height=4mm Lip=1mm of total height	TUBE22-12-4L	1.9 Guide Tube (5 Pack) ID=1.9mm OD=2.9mm Height=4mm Lip=1mm of total height	TUBE29-19-4L
1.4 Guide Tube (5 Pack) ID=1.4mm OD=2.4mm Height=4mm Lip=1mm of total height	TUBE24-14-4L	2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=4mm Lip=1mm of total height	TUBE3-2-4L
1.55 Guide Tube (5 Pack) ID=1.55mm OD=2.55mm Height=7mm	TUBE255-155-7NL	2.0 Guide Tube (5 Pack) ID=2.0mm OD=3.0mm Height=5mm Lip=1mm of total height	Tube20
1.6 Guide Tube (5 Pack) ID=1.6mm OD=2.6mm Height=4mm Lip=1mm of total height	TUBE26-16-4L	2.1 Guide Tube (5 Pack) ID=2.1mm OD=3.1mm Height=4mm Lip=1mm of total height	TUBE31-21-4L
1.7 Guide Tube (5 Pack) ID=1.7mm OD=2.7mm Height=4mm Lip=1mm of total height	TUBE27-17-4L	2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-21-5L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Size	Part #
2.1 Guide Tube ID=2.1mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-21-5L
2.1 Guide Tube ID=2.1mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-21-5L
2.11 Guide Tube (5 Pack) ID=2.11mm OD=3.11mm Height=5mm Lip=1mm of total height	Tube211
2.2 Guide Tube (5 Pack) ID=2.2mm OD=3.2mm Height=5mm Lip=1mm of total height	Tube22
2.3 Guide Tube ID=2.3mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-23-5L
2.3 Guide Tube ID=2.3mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-23-5L
2.3 Guide Tube ID=2.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-23-5L
2.31 Guide Tube (5 Pack) ID=2.31mm OD=3.31mm Height=5mm Lip=1mm of total height	Tube231

Size	Part #
2.4 Guide Tube (5 Pack) ID=2.4mm OD=3.4mm Height=4mm Lip=1mm of total height	TUBE34-24-4L
2.4 Guide Tube ID=2.4mm OD=4.0mm Height=5mm Lip High=1mm of total height Lip Wide=.5mm	TUBE4-24-5L
2.4 Guide Tube ID=2.4mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE3-24-5L
2.4 Guide Tube ID=2.4mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-24-5L
2.1 Guide Tube ID=2.1mm OD=3.0mm Height=5mm Lip=1mm of total height	TUBE4-26-5L
2.72 Guide Tube ID=2.72mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-272
2.8 Guide Tube (5 Pack) ID=2.8mm OD=3.8mm Height=5mm Lip=1mm of total height	TUBE280
2.8 Guide Tube ID=2.8mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-28-5L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Size	Part #
2.8 Guide Tube ID=2.8mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-28-5L
2.9 Guide Tube ID=2.9mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-29-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L
2.9 Guide Tube ID=2.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-29-5L
2.92 Guide Tube (5 Pack) ID=2.92mm OD=3.92mm Height=5mm Lip=1mm of total height	TUBE292
3.0 Guide Tube ID=3.0mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-3-5L
3.15 Guide Tube ID=3.15mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-315
3.3 Guide Tube ID=3.3mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-33-5L

Size	Part #
3.3 Guide Tube ID=3.3mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-33-5L
3.5 Guide Tube ID=3.5mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-35-5L
3.5 Guide Tube ID=3.5mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-35-5L
3.6 Guide Tube ID=3.6mm OD=4.0mm Height=5mm Lip=1mm of total height	TUBE4-36-5L
3.68 Guide Tube (5 Pack) 3.5mm Short/Long, Bio Trilobe ID=3.68mm OD=4.68mm Height=5mm Lip=1mm of total height	TUBE368
3.75 Guide Tube ID=3.75mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-375
3.79 Guide Tube ID=3.79mm OD=4.4mm Height=5mm Lip=1mm of total height	TUBE44-379
3.9 Guide Tube ID=3.9mm OD=5.0mm Height=5mm Lip=1mm of total height	TUBE5-39-5L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Size	Part #
4.1 Guide Tube ID=4.1mm OD=5.0mm Height=4mm Lip=1mm of total height	TUBE5-41-4L
4.1 Guide Tube (5 Pack) ID=4.1mm OD=5.1mm Height=3.5mm Lip=1mm of total height	TUBE51-41-35L
4.15 Guide Tube (5 Pack) ID=4.15mm OD=5.15mm Height=5mm Lip=1mm of total height	TUBE415
4.2 Guide Tube ID=4.2mm OD=5.0mm Height=4.5mm Lip=1mm of total height	TUBE5-42-45L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=4mm Lip=1mm of total height	TUBE525-425-4L
4.25 Guide Tube (5 Pack) ID=4.25mm OD=5.25mm Height=5mm Lip=1mm of total height	TUBE425
4.25 Guide Tube ID=4.25mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-425
4.35 Guide Tube ID=4.35mm OD=5.27mm Height=6mm Lip=1mm of total height	TUBE527-435-6L

Size	Part #
4.45 Guide Tube (5 Pack) 4.3mm Short / Long, Bio Trilobe ID=4.45mm OD=5.45mm	TUBE445
4.52 Guide Tube ID=4.52mm OD=5.57mm Height=4mm Lip=1mm of total height	TUBE557-452-4L
4.68 Guide Tube ID=4.68mm OD=5.3mm Height=5mm Lip=1mm of total height	TUBE53-468
4.7 Guide Tube (5 Pack) ID=4.7mm OD=5.7mm Height=4mm Lip=1mm of total height	TUBE57-47-4L
5.0 Guide Tube (5 Pack) ID=5.0mm OD=6.0mm Height=3.5mm Lip=1mm of total height	TUBE6-5-35L
5.0 Guide Tube ID=5.0mm OD=6.2mm Height=4mm Lip=1mm of total height	TUBE62-5-4L
5.02 Guide Tube (5 Pack) ID=5.02mm OD=6.02mm Height=5mm Lip=1mm of total height	TUBE502
5.06 Guide Tube (5 Pack) ID=5.06mm OD=6.02mm Height=3.5mm Lip=1mm of total height	TUBE6-506-35L

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



Size	Part #
5.08 Guide Tube (5 Pack) ID=5.08mm OD=6.08mm Height=5mm Lip=1mm of total height	TUBE508
5.1 Guide Tube ID=5.1mm OD=6.1mm Height=6mm Lip=1mm of total height	TUBE61-51-6L
5.15 Guide Tube (5 Pack) ID=5.15mm OD=6.15mm Height=4mm Lip=1mm of total height	TUBE615-515-4L
5.16 Guide Tube (5 Pack) 5.0mm Short, Bio Trilobe ID=5.16mm OD=6.16mm Height=5mm Lip=1mm of total height	TUBE516
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=4mm Lip=1mm of total height	TUBE625-525-4L
5.25 Guide Tube (5 Pack) ID=5.25mm OD=6.25mm Height=5mm Lip=1mm of total height	TUBE525
5.31 Guide Tube ID=5.31mm OD=6.2mm Height=4.5mm Lip=1mm of total height	TUBE62-531-45L
5.35 Guide Tube (5 Pack) ID=5.35mm OD=6.35mm Height=4mm Lip=1mm of total height	TUBE635-535-4L

Size	Part #
5.5 Guide Tube ID=5.5mm OD=6.5mm Height=4mm Lip=1mm of total height	TUBE65-55-4L
5.7 Guide Tube (5 Pack) ID=5.7mm OD=6.7mm Height=4mm Lip=1mm of total height	TUBE67-57-4L
6.13 Guide Tube ID=6.13mm OD=6.95mm Height=3.2mm Lip=1mm of total height	TUBE695-613-32L
6.2 Guide Tube ID=6.2mm OD=6.9mm Height=4mm Lip=1mm of total height	TUBE69-62-4L
6.2 Guide Tube (5 Pack) ID=6.2mm OD=7.2mm Height=3.5mm Lip=1mm of total height	TUBE72-62-35L
6.26 Guide Tube (5 Pack) ID=6.26mm OD=7.21mm Height=3.5mm Lip=1mm of total height	TUBE72-626-35L
6.3 Guide Tube ID=6.3mm OD=7.3mm Height=6mm Lip=1mm of total height	TUBE73-63-6L
6.72 Guide Tube ID=6.72mm OD=7.17mm Height=5mm Lip=1mm of total height	TUBE717-672

All prices are available in the product section of our website. Volume discounts are available for some items.

Blue Sky Bio is not affiliated with Nobel Biocare.

Blue Sky Bio's products are not manufactured, endorsed, or warranted by Nobel.



## Terms and Conditions

### Blue Sky Bio, LLC - TERMS AND CONDITIONS – Important User Information

Blue Sky Bio offers high quality dental implant products at much lower prices, compared to its competitors. To help keep our overhead low and prices value-minded, we only accept returns or exchanges under *very limited circumstances*.

#### RETURNS & EXCHANGES

Blue Sky Bio will accept returns or exchanges within 30 days of original sales receipt. Products returned or exchanged between 31 and 120 days of original sales receipt are subject to a 15% restocking fee. No products can be returned or exchanged after 120 days from date of original sales receipt. Prior to return or exchange please contact us for a Return Authorization Number (RA#). All returned or exchanged products are subject to inspection. Products that have been opened, damaged or marked upon are not eligible for return or exchange. When your products arrive, inspect the packaging and the contents for damage and defects. You must notify us within 7 days from date of original sales receipt if, product packaging is damaged, illegible, or has a broken safety seal; or if product itself appears to be damaged or does not function as its intended use. In these cases, Blue Sky Bio will pay return shipping charges, and will pay shipping charges of replacement products. Blue Sky Bio does not offer a credit or reimburse for implants that do not integrate.

#### ORDER PLACEMENT

Blue Sky Bio accepts payment for purchases with Visa, MasterCard, and American Express, through its secure website or in limited circumstances, with our fax order form. Prices for products, sales tax, shipping, or any other charges will be listed on your sales receipt. Full payment is due before orders are shipped. Sales are final after shipment has left our distribution center. Blue Sky Bio reserves the right to hold orders until credit card information can be confirmed. You may be contacted by a Blue Sky Bio representative regarding your purchase.

#### SHIPPING

Orders placed within the United States are shipped via UPS. Orders placed outside the United States are shipped via Fedex. Customers are ultimately responsible to ensure that carriers will deliver at the requested time in their areas. Orders placed before 1:30pm CST will ship the same day. Orders placed after 1:30pm CST may not ship until the next business day. Blue Sky Bio recognizes holidays throughout the year. Watch our home page for notices regarding company closings due to holidays and shipping delays due to weather.

#### INTERNATIONAL BROKER / TAX TERMS

Blue Sky Bio utilizes Fedex Express (FE) to ship and broker all orders from the United States to your door. As broker, FE pays all taxes, tariffs, and VAT due on shipments. FE will bill you separately for these charges after you receive your order, including their broker fee. Taxes and broker fees are not part of your shipping charge collected at time of order. By placing your order you authorize Blue Sky Bio to collect any taxes not paid, as soon as reported by FE of tax delinquency, via your credit card on file. You will be responsible for any additional charges or late fees associated with not paying responsible taxes, tariffs, and VAT promptly.

#### CONTACTING BLUE SKY BIO

Contacting us is very easy and convenient. Call us toll-free 9am – 5pm CST at, **888-446-6724**, or email us anytime at, [info@blueskybio.com](mailto:info@blueskybio.com).

#### IMPORTANT PRODUCT NOTICES FOR BLUE SKY BIO CUSTOMERS

BIO | Max is NobelActive® compatible.

BIO | Max implants are not compatible with Nobel Biocare Gold Adapt non-engaging UCLA abutments.

Blue Sky Bio is not affiliated with Nobel Biocare®. NobelActive® is a registered trademark of Nobel Biocare.

