

MOBILEX[®]

Lightweight & Handheld



FDA APPROVED K180561

The **AFFORDABLE** Portable Solution Safe, Effective, and Easy to Use

The MobileX produces significantly less radiation than wall mount units. Ideal for instances such as Pedo, Endo or Oral Surgery or for off site purposes when the operator must be in the room or near the patient during exposure.

Compared To National Brands, MobileX Costs Thousands Less!

The X-Ray That Goes From Op To Op, and Replaces Bulky & Expensive Wall Mounted Units!

At only 4.6 lbs. with extended cone & shield, MobileX is lighter than the competitors.

VALUE, QUALITY & SAFETY

RECHARGEABLE LI BATTERY	Take over 150 shots on one battery charge
TOSHIBA X-RAY TUBE	The best quality X-Ray tube available. Same tube as in many wall mount machines. Compare durability to lower quality tubes made in China.
OIL COOLED VS AIR COOLED	Just like your wall mount units, the MobileX features an oil-cooled tube. Oil cooling provides faster disbursement providing for faster recovery and longer lasting tubes.
FDA APPROVED	Meets all FDA standards for hand-held use. K# 180561
ATTACHED CONE & SHIELD	Mandated by most states. The attached, 20cm cone and lead infused shield protects both the operator and the patient from backscattered X-rays.
DOUBLE LEAD SHIELDED	Housing Prevents X-ray leakage from the machine and tube area.

FEATURES

- Exposure time adjustment
- Tooth chart presets
- LED back-lit screen. Displays: time, mode, battery life, ready/indicator
- Main power on/off button
- Optional remote activator switch plug receptacle with 8-foot cord
- Exposure trigger switch
- Rechargeable battery with 1-year warranty
- Extended cone. Attached cone complies with federal & state regulations
- Lead-infused back-scatter shield. Protects the operator from back-scatter

KIT INCLUDES

- Attached backscatter shield & cone
- Rechargeable battery
- Aluminum carrying case with lock & key
- Adjustable neck strap



ACTUAL IMAGE

WHY HAND-HELD VS WALL MOUNT?

Save Thousands When Equipping Your Office With MobileX Hand-Held X-Rays

- Eliminate multiple x-ray units & declutter your operatories
- Simplify your design and reduce cabinet costs
- Eliminate between operator shielding (*check with your state regulatory office*)
- Mounting: tripod compatible & the ability to use on existing wall mount arm

CONTACT YOUR DEALER TODAY!

[DEALER TAG HERE]

OPTIONAL ACCESSORIES



- Remote trigger switch w/6ft cord
- Pistol grip trigger switch

DENTERPRISE
INTERNATIONAL, INC

Denterpriseintl.com