

# Pro 95



## SprintRay Pro Dental 3D Printers



# Pro 55

1836cm (180mm x 102mm)

**Build Platform Area**

619cm (105mm x 59mm)

95um + Pixel Toning

**Minimum Feature Size**

55um + Pixel Toning

Up to 3 Inches per Hour

**Print Speed**

Up to 3 Inches per Hour

SprintRay, NextDent, Dreve, Bego, DENTCA, KeyStone

**Supported Resin Brands**

SprintRay, NextDent, Dreve, Bego, DENTCA, KeyStone

Phone Setup Consultation, Free Lifetime Phone & Email Support

**Live Support**

Phone Setup Consultation, Free Lifetime Phone & Email Support

Direct USB Printing, Contextual Help, Print Monitoring

**Onboard Computer**

Direct USB Printing, Contextual Help, Print Monitoring

\$6,750

**MSRP**

\$6,750

Full-arch Surgical Guides, Waxup Models, Retainer Models,  
Study Models, Occlusal Guards, Indirect Bonding Trays

**Recommended Indications**

Indirect Restorations, Die Models, Implant Dentures,  
Full-arch Surgical Guides, Temporaries, Cast Applications

General Practitioners, Oral Surgeons

**Target Segment**

Labs, Restorative Dentists, Prosthodontists